



WAQAR INTERNATIONAL

Care For Your Health

SURGICAL
Instruments

2018
Catalogue

www.waqarintl.com

Friendly Specialist
SERVICE

TC SCISSOR



WAQAR INTERNATIONAL
Care For Your Health



Gorney Facelift Scissors
Art# WI-S-3001
Size: 19cm
Detail: Gorney Facelift Scissors Curved TC Blades Serrated



Metzenbaum Dissecting Scissors
Art# WI-S-3002
Size: 31cm
Detail: Metzenbaum Dissecting Scissors TC Curved



Metzenbaum Dissecting Scissors
Art# WI-S-3003
Size: 31cm
Detail: Metzenbaum Dissecting Scissors TC Curved

TC FORCEPS



ADSON Dressing Forceps TC
Art# WI-S-3101
Size: 12 cm, 15 cm
Detail: ADSON Dressing Forceps TC (Tungsten Carbide)



ADSON Dressing Forceps TC
Art# WI-S-3101
Size: 12 cm
Detail: ADSON Dressing Forceps TC (Tungsten Carbide)
1x2 teeth Platform cross serrated



ADSON Dressing Forceps TC
Art# WI-S-3101
Size: 12 cm, 15 cm
Detail: ADSON Dressing Forceps TC (Tungsten Carbide)
1x2 teeth Serrated jaws

BONE RONGEUR



WAQAR INTERNATIONAL
Care For Your Health



Beyer Bone Rongeur

Art# WI-2031

Size: 180mm

Detail: End-cutting rongeurs used for trimming and recontouring alveolar bone and gross tissue removal. Slightly, Curved.



Friedman Bone Rongeur

Art# WI-2032

Size: 140mm, 120mm

Detail: End-cutting rongeurs used for trimming and recontouring alveolar bone and gross tissue removal. Slightly, Curved.



Kerrison Bone Rongeur

Art# WI-2033

Size: 30mm

Detail: End-cutting rongeurs used for widening the sinus window.

CARDIOVASCULAR CLAMPS



Aortic Exclusion Clamp

Art# WI-1950

Size: 8.5 cm, 7-1/4" (18.5 cm)

Detail: The unique "S" shaped curvature of the jaws makes it contoured to fit.



Castaneda Clamp

Art# WI-1951

Size: 6" (15.0 cm)

Detail: The curved shanks, long thin jaw, calibrated jaws, horizontal serrations.



Cooley Anastomosis Clamp

Art# WI-1952

Size: 6-1/2" (16.5 cm)

Detail: The jaws calibrated at 5.0 mm intervals



Gregory Profunda Clamp

Art# WI-1953

Size: 7-3/4" (19.5 cm)

Detail: The large jaws.


SELF RETAINING RETRACTORS



WAQAR INTERNATIONAL
Care For Your Health



Henly Carotid Retractor
Art# WI-1132
Size: 6-1/2" (16.5 cm)
Detail: Modified weitlaner frame
Center blades, Stainless Steel



Parsonnet Epicardial Retractor
Art# WI-1133
Detail: 3x3 sharp prongs
Available Sizes are given below:
1-1/8" (32.0 mm)
1-1/2" (38.0 mm)
1-7/8" (48.0 mm)



Kerrison Bone Rongeur
Art# WI-2033
Detail: Direction of widest blade upon placement determines whether Right
Length: 6-1/2" (16.5 cm)
Spread: 3-3/4" (9.5 cm)
Large blade: 1-1/2" (3.8 cm) Wide x 3" (7.5 cm) Deep
Small blade: 1" (2.5 cm) Wide x 2-1/8" (5.4 cm) Deep
Right

SUCTION TUBES



Cooley Sump Suction Tube
Art# WI-2511
Size: 7.5 mm, 13" (33.0 cm)
Detail: Wide tip, entire tube removable from handle, 13" (33.0 cm) working length



Baron Suction Tubes Set of 3
Art# WI-2512
Detail:
Baron Suction Tube #3FR, Size 1mm.
Baron Suction Tube #5FR, Size 1.7mm.
Baron Suction Tube #7FR, Size 2.3mm

Dental Sinus Suction Tubes Set of 3
Art# WI-2513
Detail:
Suction Tube 2.5mm diameter. 6 inches long.
Suction Tube 3mm diameter. 6 inches long.
Suction Tube 4mm diameter. 6 inches long.

TWEEZERS



WAQAR INTERNATIONAL
Care For Your Health



Surgical Forceps Ewald

Art# WI-3201

Size: 13 cm

Detail: Surgical instruments are European made high quality products of stainless steel.



Surgical Suture Forceps

Art# KT-3202

Size: 13 cm

Detail: German quality stainless steel suture forceps.



Towel clip Jones

Art# WI-3203

Size: 9 cm

Detail: Standard quality stainless steel towel clip Jones.



Ear Forceps

Art# WI-3204

Size: 14 cm

Detail: Standard quality stainless steel towel ear Forceps.

MALLETS, TAMPS, IMPACTORS



SOLID COPPER HEAD MALLET

Art# KT-2210

Size: 3/4" (20mm), 6 1/2" (16.5cm)

Detail: Solid COPPER HEAD MALLET, Head diameter 3/4" (20mm), 6 1/2" (16.5cm)



NYLON HEAD MALLET

Art# WI-2211

Size: 25.4mm, 7.76 oz., 7 1/2" (19.0cm)

Detail: NYLON HEAD MALLET, Head diameter 25.4mm, 7.76 oz., 7 1/2" (19.0cm)

MEAD MALLET

Art# WI-2211


Size: 7 1/2" (19.0cm)

Detail: MEAD MALLET, Head diameter 26.0mm, Weight 11 oz (315 grams), Head Weight: 6 oz (170 grams),


PLASTER INSTRUMENTS



WAQAR INTERNATIONAL
Care For Your Health



Plaster Instruments
Art# WI-2403
Size: 170 mm
Detail: Jaws 17-198 with longitudinal and transversal grooves




Plaster Instruments
Art# WI-2411
Size: 1.8mm
Detail: Bending Plier For Wire Up to 1.8mm diam.



Plaster Instruments
Art# WI-2426
Size: 1.8mm
Detail: Cutter for cutting off KIRSCHNER wires, Tibia head screws, bolts, pins, up to 3 mm diam. Tibia head screws, bolts, pins, up to 6 mm diam.


WOUND RETRACTORS



Wound Retractors
Art# WI-1002
Size: 125 mm



Wound Retractors
Art# WI-1003
Size: 125 mm



Wound Retractors
Art# WI-1006
Size: 135 mm



Wound Retractors
Art# WI-1016
Size: 155 mm


WIRE CUTTING PLIER WITH T C INSERTS




WAQAR INTERNATIONAL
Care For Your Health



Wire Cutting Plier With T C Inserts
Art# WI-2901
Size: 17.5 cm
Detail: Plier | 17.5 cm



Wire Cutting Plier With T C Inserts
Art# WI-2902
Size: 17.5 cm
Detail: Plier | 17.5 cm



Wire Cutting Plier With T C Inserts
Art# WI-2903
Size: 17.5 cm
Detail: Plier | 17.5 cm

SCALPELS KNIVES AND SCALPEL HANDLES




Scalpels, Knives And Scalpel Handles
Art# WI-201
Size: 27 mm



Scalpels, Knives And Scalpel Handles
Art# WI-202
Size: 35 mm



Scalpels, Knives And Scalpel Handles
Art# 203
Size: 47 mm




Scalpels, Knives And Scalpel Handles
Art# 204
Size: 55 mm


NEEDLE HOLDERS, SUTURE INSTRUMENTS




VAQAR INTERNATIONAL
Care For Your Health



Needle Holders, Suture Instruments
Art# WI-816
Size: 150, 180, 200 cm
Detail: MAYO-HEGAR
medium fine patterns



Needle Holders, Suture Instruments
Art# WI-817
Size: 150, 180, 200 cm
Detail: MAYO-HEGAR
medium fine patterns



Needle Holders, Suture Instruments
Art# WI-818
Size: 150, 180, 200 cm
Detail: MAYO-HEGAR
medium fine patterns

SUCTION INSTRUMENTS



Suction Instruments
Art# WI-2501



Suction Instruments
Art# WI-2502



Suction Instruments
Art# WI-2504






Suction Instruments
Art# WI-2508

- 18-167 Ø 2,0 mm
- 18-168 Ø 2,5 mm
- 18-169 Ø 3,0 mm
- 18-170 Ø 4,0 mm
- 18-171 Ø 5,0 mm

BONE PUNCHES RONGEURS



WAQAR INTERNATIONAL
Care For Your Health



Bone Punches Rongeurs

Art# WI-2001

Size: 180 mm

Detail: Kerrison Jacoby - Ferris Smith
Kerrison Angled 130 Degree upwards
cutting Shafts Sizes: "180 mm"
Punches Diam

Bone Punches Rongeurs

Art# WI-2002

Size: 180 mm

Detail: Kerrison Jacoby - Ferris Smith
Kerrison Angled 130 Degree upwards
cutting Shafts Sizes: "180 mm"
Punches Diam


Bone Punches Rongeurs

Art# WI-2003

Size: 180 mm

Detail: Kerrison Jacoby - Ferris Smith
Kerrison Angled 130 Degree upwards
cutting Shafts Sizes: "180 mm"
Punches Diam

SUCTION INSTRUMENTS



Suction Instruments

Art# WI-2501



Suction Instruments

Art# WI-2502



Suction Instruments

Art# WI-2504



Suction Instruments

Art# WI-2508

- 18-167 Ø 2,0 mm
- 18-168 Ø 2,5 mm
- 18-169 Ø 3,0 mm
- 18-170 Ø 4,0 mm
- 18-171 Ø 5,0 mm



WAQAR INTERNATIONAL

Care For Your Health

Waqar International

Address: New Miana Pura, Gulshan Town,
Near Railway Line Sialkot - 51310 Pakistan.

Call us: +92 304 4300461

Email us: info@waqarintl.com