

WAQAR INTERNATIONAL

Precision Instruments Manufacturer

Premium Quality Instruments Manufacturer Of Surgical, Dental, Veterinary Instruments.



Our Range of Best Selling Instruments

SURGICAL & DENTAL

INSTRUMENTS CATALOGUE



INSTRUCTIONS FOR REPROCESSING OF REUSABLE SURGICAL INSTRUMENTS

The following instructions are for all reusable medical devices supplied by WAQAR INTERNATIONAL unless stated otherwise with the packaging of the product. These instructions are intended for use only by persons with the required specialist knowledge and training.

Warnings

- Follow instructions and warnings as issued by the manufacturers of any decontaminants, disinfectants and cleaning agents. Wherever possible avoid use of mineral acids and harsh, abrasive agents.
- No part of the process shall exceed 140 °C.
- Some sensitive materials (e.g. Aluminium) are damaged by high alkaline solutions (pH>10).
- Devices with long, narrow cannula, hinges and blind holes require particular attention during cleaning.

Note: When reprocessing medical devices, always handle with care, wearing protective clothing, gloves and eye-wear in accordance with local health & safety procedures.

Limitations on Reprocessing

- Repeated processing has minimal effect on these instruments.
- End of life normally determined by wear and damage in use.
- Any specific limitations on the number of reprocessing cycles shall be made available with instrument.

From Point of Use

Whenever possible, do not allow blood, debris or bodily fluids to dry on instruments. For best results and to prolong the life of the medical device reprocess immediately after use. If they cannot be reprocessed immediately, use an enzymatic foam spray cleaner to help prevent soil from drying.

Preparation for Decontamination

- Process all the instruments as soon as it is reasonably practical following use.
- Disassemble only where intended, without the use of tools unless specifically provided by the manufacturer. Instructions for dis assembly are available with the device, where required.

Cleaning: Automated

- Use only either CE marked or validated washer-disinfector machines and low foaming, non-ionizing cleaning agents and detergents following the manufacturer's instructions for use, warnings, concentrations and recommended cycles.
- Load instruments carefully, with any box joints and hinges open and so that any fenestration in instruments can drain.
- Place heavy instruments with care in the bottom of containers, taking care not to overload wash baskets.
- Place instruments with concave surfaces facing down to prevent pooling of water.
- Where possible, use appropriate attachments to flush inside reamers and devices with lumens or cannula.
- Ensure that soft, high purity water which is controlled for bacterial endo-toxins is used in the final rinse stage.

Note: Automated cleaning may not be suitable for all lumens and cannula in which case clean manually with a water jet gun, if available, and an appropriate brush (and stilette if provided) that reaches the depth of the feature. After manually cleaning, pass all devices through an automatic cleaning cycle to achieve disinfection.

Note: These instruments have been validated using a washer-disinfector cycle validated to include two cold rinses at <35 °C, a detergent cycle and a rinse cycle both at >50°C, a disinfection cycle operating at a temperature of between 80 °C and 90 °C for a minimum holding time of 1 minute (actual holding time in excess of 2 minutes 50 seconds) and a 20 minute drying cycle. The detergent used was low foaming, non-ionizing spray wash detergent cleaner (max 12pH) and the rinse aid a neutral pH low foaming, non-ionic surfactant with isopropyl alcohol.

Cleaning: Manual

- Manual cleaning is not advised if an automatic washer-disinfector is available. If the equipment is not available, use the following process:-
1. Use a double sink system (wash/rinse) dedicated for instrument cleaning (not used for hand washing). Ensure that the water temperature does not exceed 35 of the re-processor to ensure that the reprocessing as actually performed using equipment, materials and personnel in the reprocessing facility achieve the desired result. This requires validation and routine monitoring of the process. Likewise any deviation by the re-processor from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences.
 2. In the first sink, keeping the instrument submerged, with an autoclavable brush, apply CE marked cleaning solution to all surfaces until all soil has been removed. Pay particular attention to serrations, teeth, ratchets and hinges, always brushing away from the body and avoiding splashing. Ensure rongeurs and hinged instruments are thoroughly cleaned in both open and closed positions.
 3. In the second sink, rinse instruments thoroughly with soft, high purity water which is controlled for bacterial endotoxins, so that the water reaches all parts of the instrument, then carefully hand dry or use a drying cabinet.

Note: Manual Cleaning is NOT a disinfection process: When manual cleaning is used it may not be possible to disinfect the device prior to further handling.

Cleaning Inspection

- After cleaning visually inspect, visually inspect all surfaces, cannulations, ratchets, joints, holes and lumens for complete removal of soil and fluids. If ANY soil or fluid is still visible, return the instrument for repeat decontamination.

Maintenance

- Apply surgical grade lubricant to hinges, joints and moving parts as per the lubricant manufacturer's instructions.

Continued on next page

Inspection and Functional Testing

- Visually inspect and check: - all instruments for damage and wear; cutting edges are free from nicks and present a continuous edge; jaws and teeth align correctly; all articulated instruments have a smooth movement without excess play; locking mechanisms (such as ratchets) fasten securely and close easily; long, slender instruments are not distorted; any component parts fit and assemble correctly with mating components.
- Remove for repair and replacement any blunt, worn out, flaking, fractured or damaged instruments.

Note: If an instrument is returned to the manufacturer / supplier, the instrument must be decontaminated and sterilized and be accompanied with the relevant documented evidence.

Packaging

- All instruments to be packed following local protocol in accordance with ISO 11607-1 or AAM/CSR technique.

Sterilization (USA)

- Autoclave should comply with the requirements of, and be validated and maintained in accordance with ANSI/AAMI ST79.
- Pre-Vacuum moist heat sterilization operating at temperature 132 °C bar for a minimum holding time of 4 minutes always following the instructions of the machine manufacturer.
- When sterilizing multiple instruments in one autoclave cycle, ensure that the sterilizer manufacturers stated maximum load is not exceeded.
- Ensure instruments are dry before sterilization.

Sterilization (Outside USA)

- Autoclave should comply with the requirements of, and be validated and maintained in accordance with EN 285, EN 13060, EN ISO 17665.
- Either CE marked or validated vacuum autoclave operating at 134-137 °C bar for a minimum holding time of 3 minutes always following the instructions of the machine manufacturer.
- When sterilizing multiple instruments in one autoclave cycle, ensure that the sterilizer manufacturers stated maximum load is not exceeded.

Storage

- Ensure instruments are dry before storage, and stored in dry, clean conditions at an ambient room temperature.

Additional Information

- Other forms of cleaning (i.e. ultrasonic) and sterilization (i.e. Low temperature steam and Formaldehyde, Ethylene Oxide and Gas Plasma) are available. However, always follow the instructions for use as issued by the manufacturer and always consult with them if in any doubt over the suitability of any process used.
- Cleaning and sterilizing guidelines are available in HTM 2030 and HTM 2010. Contact: The NHS Estates Stationers Office Publications Cent for details at www.tsonline.gov.uk. For further information contact: NHS Estates Information Center, Department of Health, 1 Trevelyan Square, Boar Lane, Leeds, LS1 6AE, UK or visit www.nhsestates.gov.uk.

Note: It is the responsibility of the re-processor to ensure that the reprocessing as actually performed using equipment, materials and personnel in the reprocessing facility achieve the desired result. This requires validation and routine monitoring of the process. Likewise any deviation by the re-processor from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences.





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Cesta para Residuos
Satinado, Satin

Catalog # _____

WI-010-03



15 L

Cheattle Sterilizing Forceps
Satin

Catalog # _____

WI-010-04

Size _____

27 cm – 10 3/4



Clip de Esterilización de Mayo
Satinado, Satin

Catalog # _____

WI-010-34 14 cm – 5 1/2 , 14 cm – 5 1/2

Size _____



Plain End



Ball End

Deponette Storing Case
Mirror

Catalog # _____

WI-012-03



Perforated Bottom with Paper Filter

Harrison Sterilizing Bowl Forceps
Satin

Catalog # _____

WI-020-13




Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris
 Manici per bisturi

Martin Waste Forceps
 Satin



Catalog #	Size
WI-020-14	58 cm – 22 3/4

Mayo Safety Pin
 Satin



Catalog #	Size
WI-021-13	14 cm – 5 1/2

Pinza para Residuos de Sauerbruch
 Satinado, Satin



Catalog #	Size
WI-021-14	55 cm – 21 1/2 , 55 cm – 21 1/2

Rogge Sterilizing Forceps
 Satin



Catalog #	Size
WI-030-07	30.5 cm – 12

Round Needle Case
 Satin



Complete with Perforated Interior box



Perforated Lid and Bottom

Catalog #	Size
WI-030-17	Ø 65 x 15 mm



Scalpel Handles
Skalpellgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Sauerbruch Needle Case
Satin

Catalog # _____ Size _____
WI-020-14 210 x 125 x 17 mm



Sauerbruch Waste Forceps
Satin

Catalog # _____ Size _____
WI-021-13 55 cm – 21 1/2



Schimmelbusch Sterilization Drum
Mirror

Catalog # _____
WI-021-14



Silicone Handles

Catalog # _____
WI-030-07




Sterilization Boxes
Mirror

Catalog # _____
WI-030-17



Perforated Lid and Bottom



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi


Sterilization Cassettes
 Mirror

Catalog # _____
 WI-020-14




Surticon Dividers
 Mirror

Catalog # _____
 WI-021-13



Surticon Drape Retainers
 Mirror

Catalog # _____
 WI-021-14




Surticon Identification Labels

Catalog # _____
 WI-030-07



Surticon Instrument Basket Stamp

Catalog # _____
 WI-030-17



20 mm


14 mm



Scalpel Handles
Skalpelli-griffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi


Surticon Instrument Baskets
Mirror

Catalog # _____
WI-020-14



Surticon Instrument Racks
Mirror

Catalog # _____
WI-021-13



Surticon Label Holder
Mirror

Catalog # _____ Overall Size
WI-021-14 50 x 35 x 13



Surticon Security Seal

Catalog # _____
WI-030-07





With Indicator (Pack of 100), Blue



With Indicator (Pack of 100), Red

Surticon Silicone Mats


Catalog # _____
WI-030-17



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
ici per bisturi

Surticon Single Use Paper Filters

Catalog # _____
WI-020-14



Surticon Teflon Filters

Catalog # _____
WI-021-13



Surticon Wire Mesh Baskets
Mirror

Catalog # _____
WI-021-14



Textile Filters

Catalog # _____
WI-030-07




Package of 5 Pairs

Waste Basket
Satin

Catalog # _____ Overall Size
WI-030-17 325 x Ø 305 mm



15 L




Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Weinstein Instrument Holder
Satin

Catalog #

WI-020-14



Adams Napkin Holder
Satin

Catalog #
WI-020-14




Consisting of Chain and Two Crocodile Forceps

Scalpel Handles
Skalpgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Agarwal Karate Phaco Chopper
Satin

Catalog # Overall Length
WI-021-13 11.5 cm – 4 1/2



Akahoshi Nucleus Sustainer
Satin

Catalog # Overall Length
WI-021-14 11.5 cm – 4 1/2



Alfonso Cutting Platform
Satin

Catalog # Overall Length
WI-030-07 12.5 cm – 5



Alfonso Nucleus Trisector
Satin

Catalog # Overall Size
WI-030-17 10.5 cm – 4 1/4





Scalpel Handles
Skalpelligriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Amalgam Carrier
Satin

Catalog # _____

WI-020-14



Amalgam Well
Satin

Catalog # _____

WI-021-13



Angle Former Fig. 30/31
Satin

Catalog # _____

WI-021-14



Angle Former Fig. 32/33
Satin

Catalog # _____

WI-030-07



Angle Former Fig. 34/35
Satin

Catalog # _____

WI-030-17



Scalpel Handles
 Skalpelligriffe
 Mangos de bisturios
 Manches pour bistouris
 Manici ner histuri

Anis Nucleus Manipulator
 Satin

Catalog #	Overall Length
WI-020-14	12 cm – 4 3/4



Asa de New Orleans
 Satin

Catalog #	Overall Length
WI-021-13	13 cm – 5 , 13 cm – 5



Asa de Wilder
 Satin

Catalog #	Overall Length
WI-021-14	13.5 cm – 5 1/4 , 13.5 cm – 5 1/4



Aspiration Uterine Curette and Coupling Handle Piece
 Satin

Catalog #	Overall Length
WI-030-07	29 cm – 11 1/2



Babinsky Percussion Hammer
 Satin

Catalog #	Overall Size
WI-030-17	22 cm – 8 3/4





Scalpel Handles
Skalpellgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Baldwin Plastic Filling Instrument Fig. 49
Satin

Catalog #

WI-020-14



Balshi Gingival Cord Packer
Satin

Catalog #

WI-021-13



Band Seater and Pusher
Satin

Catalog #

WI-021-14



Bayonet Shaped Mirror Handle
Satin

Catalog #

WI-030-07

Overall Length

15 cm – 6



Bayonet Shaped

Bechert Nucleus Rotator
Satin

Catalog #

WI-030-17

Overall Length

11.5 cm – 4 1/2



Blunt Forked Tip
Angled at 7 mm



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

BEQ1 Quetin Plastic Filling Instrument
 Satin



Catalog #

WI-020-14

Boghosian Plastic Filling Instrument Fig. Ab1
 Satin



Catalog #

WI-021-13

Boghosian Plastic Filling Instrument Fig. Ab1
 Satin



Catalog #

WI-021-14

Boghosian Plastic Filling Instrument Fig. Ab2
 Satin



Catalog #

WI-030-07

Bohler Extension Bow
 Satin



Catalog #

WI-030-17



Scalpel Handles
Skalpelligriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Bone Graft Syringe
Satin



Catalog #

WI-020-14

Bone Packer
Satin



Catalog #

Working End Size

WI-021-13

5.0 mm / 7.0 mm

Brand Tendon Tunnelling Forcep
Satin



Catalog #

Working End Details

WI-021-14

Serrated Jaws with 1 x 2 Teeth

Braun Uterine Depressor
Satin



Catalog #

Overall Length

WI-030-07

27 cm – 10 3/4

Brierley Nucleus Splitter
Satin



Catalog #

Overall Length

WI-030-17

12 cm – 4 3/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturias
 Manches pour bistouris
 Manici per bisturi

Brierley Nucleus Splitter
 Satin



Catalog #	Overall Length
WI-020-14	12 cm – 4 3/4

Brompton Tube Milking Rolle
 Blue Anodized



Catalog #	Overall Length
WI-021-13	

Bruenings Tonsil Snare
 Satin



Catalog #	Overall Length
WI-021-14	29 cm – 11 1/2

Bucher-Kansas Nucleus Trisector
 Satin



Catalog #	Overall Length
WI-030-07	10.5 cm – 4 1/4

Buck Percussion Reflex Hammer
 Satin



Catalog #	Overall Length
WI-030-17	18.5 cm – 7 1/4



Scalpel Handles
Skalpgriffe
Manches pour bistouris
Mangos de bisturíes
Manici per bisturi



Bunge Evisceration Spoon
Satin

Catalog #	Overall Length
WI-020-14	14 cm – 5 1/2



Bunnell Tendon Passer
Satin

Catalog #	Overall Length
WI-021-13	23 cm – 9



Bunnell Tendon Stripper
Satin

Catalog #	Overall Length
WI-021-14	15 cm – 6



Burnisher Fig. 21
Satin

Catalog #
WI-030-07



Calcium Hydroxide Placement
Satin

Catalog #
WI-030-17



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manici per bisturi

Carter Sphere Introducer
 Satin

Catalog #	Overall Length
WI-020-14	12.5 cm – 5



Caspar Distraction Screw
 Satin

Catalog #	Overall Length
WI-021-13	14 mm, 16 mm, 18 mm



14 mm, 16 mm, 18 mm

Caspar Screwdriver
 Satin

Catalog #	Overall Length
WI-021-14	20 cm – 8



Cement Mixing Slab

Catalog #	Overall Length
WI-030-07	11.5 cm – 4 1/2



Chang Micro Finger Chopper
 Satin

Catalog #	Overall Length
WI-030-17	11.5 cm – 4 1/2





Scalpel Handles
Skalpelligriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Chang Phaco Chopper
Satin

Catalog #	Overall Length
WI-020-14	13 cm – 5



Choper de Hwang
Satinado, Satin

Catalog #	Overall Length
WI-021-13	13.5 cm – 5 1/4 , 13.5 cm – 5 1/4



Clayman Nucleus Rotator
Satin

Catalog #	Overall Length
WI-021-14	12 cm – 4 3/4




Cone Socket Mirror Handle
Satin

Catalog #	Overall Length
WI-030-07	



Contact Forming Instrument
Satin

Catalog #	Overall Length
WI-030-17	



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Contouring Instrument
 Satin

Catalog #
 WI-020-14



Cook Waite Dental Aspirating Syringe
 Satin

Catalog #
 WI-021-13



Cook Waite Dental Aspirating Syringe
 Satin

Catalog #
 WI-021-14



Cucharilla de Cullen
 Satinado, Satin

Catalog #	Overall Length
WI-030-07	13 cm – 5 , 13 cm – 5



Cucharilla de Paton
 Satinado, Satin

Catalog #	Overall Length
WI-030-17	14 cm – 5 1/2 , 14 cm – 5 1/2





Scalpel Handles
Skalpellgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Culler Lens Spoon
Satin



Catalog #	Overall Length
WI-020-14	13 cm – 5

Dastoor Pupil and Iris Repositor
Satin




Catalog #	Overall Length
WI-021-13	13 cm – 5

Dental Aspirating Syringe
Satin



Catalog #	Overall Length
WI-021-14	

Dodick-Kammann Nucleus Chopper
Satin



Catalog #	Overall Length
WI-030-07	12 cm – 4 3/4

Drysdale Nucleus Manipulator
Satin



Catalog #	Overall Length
WI-030-17	11 cm – 4 1/4



Ear Syringe
Satin


Catalog # _____
WI-020-14



Scalpel Handles
Skalpellgriffe
Mangos de bisturías
istouris
i

Ellis Foreign Body Spud
Satin


Catalog # _____ Overall Length _____
WI-021-13 11 cm – 4 1/4



Lightly Curved and Rounded Tip

Escaf Nucleus Divider
Satin

Catalog # _____ Overall Length _____
WI-021-14 12 cm – 4 3/4



Wedge Shaped

Escaf Nucleus Manipulator
Satin

Catalog # _____ Overall Length _____
WI-030-07 12 cm – 4 3/4



Angled 45°

Escaf Nucleus Sustainer
Satin

Catalog # _____ Overall Length _____
WI-030-17 12 cm – 4 3/4



Smooth Ball Tip
Texture interior
Surface



Scalpel Handles
Skalpellsgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Gauze Applicator
Satin

Catalog #

WI-020-14



Hartmann (French) Tuning Fork
Satin

Catalog #


WI-021-13



Hartmann Tuning Fork
Satin

Catalog #

WI-021-14



Ivory Amalgam Carrier
Satin

Catalog #


WI-030-07



Koenig Tuning Fork
Satin

Catalog #

WI-030-17



Scalpel Handles
Skalpelligriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Lucae Tuning Fork
Satin

Catalog #
WI-020-14



McCoy Laryngoscope
Satin

Catalog #
WI-021-13



McGivney Hemorrhoidal Ligator
Satin

Catalog #
WI-021-14



McGivney Suction Hemorrhoidal Ligator
Satin

Catalog # Overall Length
WI-030-07 29 cm – 11 1/2



McIntosh Laryngoscope
Satin

Catalog #
WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Membrane Placement Instrument
Satin

Catalog #

WI-020-14



Mershon Band Pusher
Satin

Catalog #

WI-021-13



Michel Laryngeal Mirror
Satin

Catalog #

WI-021-14

Overall Length

25 cm – 9 3/4



With Handle , Moveable Mirror

Microchoper de Chang
Satin

Catalog #

WI-030-07

Overall Length

11.5 cm – 4 1/2 , 11.5 cm – 4 1/2



Miller Laryngoscope
Satin

Catalog #

WI-030-17



WI-030-17 Complete Set

With Fiber Optic Complete Set

55 mm
Fig. 0
For Babies

80 mm
Fig. 1
For Children

130 mm
Fig. 2
For Adolescents

170 mm
Fig. 3
For Women

180 mm
Fig. 4
For Men

Battery Handle, Chrome Plated




Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Minardi Phaco Chopper
 Satin



Catalog #	Overall Length
WI-020-14	11.5 cm – 4 1/2

Mirror
 Satin



Catalog #	
WI-021-13	

Ø 6.0 mm, Ø 8.0 mm, Ø 10.0 mm, Ø 12.0 mm

Neuhann Nucleus Divider
 Satin



Catalog #	Overall Length
WI-021-14	12 cm – 4 3/4

Neuro Tuning Fork
 Satin



Catalog #	
WI-030-07	

With Hammer, Pinwheel, Brush and Needle

New Orleans Lens Loop
 Satin



Catalog #	Overall Length
WI-030-17	13 cm – 5



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Otis Urethrotome
Satin

Catalog #

WI-020-14



P.K Thomas Wax Instruments Fig. 2
Satin

Catalog #

WI-021-13



Fig. 2

P.K Thomas Wax Instruments Fig. 3
Satin

Catalog #

WI-021-14



Fig. 3

P.K Thomas Wax Instruments Fig. 4
Satin

Catalog #

WI-030-07



Fig. 4

P.K Thomas Wax Instruments Fig. 5
Satin

Catalog #

WI-030-17



Fig. 5



Scalpel Handles
 Skalpelligriffe
 Mangos de bisturías
 Manches pour bistouris
 Manici per bisturi

Paton Spatula & Spoon
 Satin



Catalog #	Overall Length
WI-020-14	14 cm – 5 1/2

Percy Amputation Shield
 Satin



Detachable

Catalog #
WI-021-13

Philips Screw Driver
 Satin



Catalog #	Overall Length
WI-021-14	25 cm – 10

Pisacano Nucleus Rotator
 Satin



Catalog #	Overall Length
WI-030-07	11.5 cm – 4 1/2

Plastic Filling Instrument Fig. G
 Satin



Catalog #
WI-030-17



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Plastic Filling Instrument Fig. W1
Satin



Catalog # _____

WI-020-14

Plastic Filling Instrument Fig. W3
Satin



Catalog # _____

WI-021-13

Rabiner Percussion Hammer
Satin



Catalog # _____

Overall Length

WI-021-14

23 cm – 9

Remington Hoe
Satin



Catalog # _____


WI-030-07

Roach Wax Carver
Satin



Catalog # _____

WI-030-17



Scalpel Handles
Skalpellgriffe

Rosen Phaco Splitter
Satin

Catalog #	Overall Length
WI-020-14	11.5 cm – 4 1/2



Rotador de Clayman
Satinado, Satin

Catalog #	Overall Length
WI-021-13	12 cm – 4 3/4 , 12 cm – 4 3/4



Rotador de Hunkeler
Satin

Catalog #	Overall Length
WI-021-14	12 cm – 4 3/4 , 12 cm – 4 3/4



Rotador de Jaffe-Bechert
Satin, Satin

Catalog #	Overall Length
WI-030-07	11.5cm – 4 1/2 , 11.5cm – 4 1/2



Rotador de Pisacano
Satinado, Satin

Catalog #	Overall Length
WI-030-17	11.5 cm – 4 1/2 , 11.5 cm – 4





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Rowen Combination Lasik Flap Protector and Spatula
Satin




Catalog #	Overall Length
WI-020-14	13 cm – 5

Rubber Bowl




Catalog #
WI-021-13

Rubber Dam Frame
Satin



Catalog #	Overall Length
WI-021-14	23 cm – 9

Rudd Hemorrhoidal Ligator
Satin



Catalog #	Overall Length
WI-030-07	31.5 cm – 12 1/2

Rumel-Belmont Tourniquet
Satin



Catalog #
WI-030-17

CV-794-27 27 cm, Ø 3.2 mm
CV-796-30 30.5 cm, Ø 6.4 mm



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris

Rydel-Seiffert Tuning Fork
 Satin

Catalog # _____
 WI-020-14



Schamberg Comedone Extractor
 Satin

Catalog # _____ Overall Length _____
 WI-021-13 14.5 cm – 5 3/4



Screw Driver
 Satin

Catalog # _____ Overall Length _____
 WI-021-14 25 cm – 10



Shaw Spatula Fig. 7
 Satin

Catalog # _____
 WI-030-07



Shaw Wax Instrument Fig. 1
 Satin

Catalog # _____
 WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Shaw Wax Instrument Fig. 2
Satin

Catalog #
WI-020-14




Siepsner Endocapsular Controller
Satin

Catalog # Overall Length
WI-021-13 11.5 cm – 4 1/2



Sims Uterine Depressor
Satin

Catalog # Overall Length
WI-021-14 26 cm – 10 1/4




Sloane LASEK Epi Peeler
Satin

Catalog # Overall Length
WI-030-07 11.5 cm – 4 1/2



Sloane LASEK Flap Repositor
Satin

Catalog # Overall Length
WI-030-17 12 cm – 4 3/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris


Sloane LASEK Micro Hoe
 Satin

Catalog #	Overall Length
WI-020-14	12 cm – 4 3/4



Sloane LASEK Micro Hoe & Epi Peeler
 Satin

Catalog #	Overall Length
WI-021-13	13.5 cm – 5 1/4



Slotted Ring Pattern Tendon Stripper
 Satin

Catalog #	Overall Length
WI-021-14	28.5 cm – 11 1/4



Sluder-Ballenger Tonsillectome
 Satin

Catalog #	Overall Length
WI-030-07	



Snare Wire
 Satin

Catalog #	Overall Length
WI-030-17	10.5 cm – 4





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Snare Wire
Satin

Catalog # _____
WI-020-14



Pack of 100

Snare Wire
Satin

Catalog # _____ Overall Length _____
WI-021-13 100 cm – 39 1/2



Ø 0.3 mm, Ø 0.4 mm

Sostenedor de Núcleo de Escaf
Satinado, Satin

Catalog # _____ Overall Length _____
WI-021-14 12 cm – 4 3/4 , 12 cm – 4 3/4



Suction Hemorrhoidal Ligator
Satin

Catalog # _____
WI-030-07



SI-455-01 Aligned 33° Up, With 10 mm Suction Head and Charging Cone
SI-457-01 Aligned 33° Down, With 11 mm Suction Head and Charging Cone

T-Bar Handle Bone Tap
Satin

Catalog # _____
WI-030-17



Scalpel Handles
 Skalpgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Taylor Percussion Hammer
 Satin



Catalog #	Overall Length
WI-020-14	20 cm – 8

Tendon Interlacing Forceps
 Satin



Catalog #	Overall Length
WI-021-13	13.5 cm – 5 1/4

Tendon Pulling Instrument
 Satin



Catalog #	Overall Length
WI-021-14	33.5 cm – 13 1/4

Tendon Stripper
 Satin



Catalog #	Overall Length
WI-030-07	23 cm – 9

Tennant Nucleus Rotator
 Satin



Catalog #	Overall Length
WI-030-17	12 cm – 4 3/4



Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Torniquete de Feldheim
Satinado, Satin

Catalog #	Overall Length
WI-020-14	10 cm – 4 , 10 cm – 4



Trico Plastic Filling Instrument Fig. 1/2
Satin

Catalog #	Overall Length
WI-021-13	100 cm – 39 1/2



Trisector de Kansas
Satinado, Satin

Catalog #	Overall Length
WI-021-14	11 cm – 4 1/4 , 11 cm – 4 1/4



Tufts Plastic Filling Instrument Fig. 3
Satin

Catalog #	Overall Length
WI-030-07	



Tufts Plastic Filling Instrument Fig. 4F
Satin

Catalog #	Overall Length
WI-030-17	



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi




Tuning Fork
 Satin

Catalog # _____
 WI-020-14

Tunneling Forcep
 Satin

Catalog # _____ Overall Length _____
 WI-021-13 21 cm – 8 1/4



Tyding Tonsil Snare
 Satin


Catalog # _____
 WI-021-14



With 1 Straight Tip

Unna Comedone Extractor
 Satin

Catalog # _____ Overall Length _____
 WI-030-07 14.5 cm – 5 3/4



Vehe Wax Carver
 Satin

Catalog # _____
 WI-030-17






Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Vollmar Ring Stripper
Satinado, Satin

Catalog #	Overall Length
WI-020-14	50 cm – 19 3/4



Wax Spatula
Mirror

Catalog #	
WI-021-13	



Wax Spatula Fig. 31
Mirror

Catalog #	
WI-021-14	



Wax Spatula Fig. 7
Mirror

Catalog #	
WI-030-07	



Wells Enucleation Spoon
Satin

Catalog #	Overall Length
WI-030-17	15 cm – 6




Scalpel Handles
 Skalpgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Wilder Lens Loop
 Satin



Catalog #	Overall Length
WI-020-14	13.5 cm – 5 1/4

Wire Loop for Tonsil Snare
 Satin




Catalog #
WI-021-13

Wire Strainer
 Satin



Catalog #	Overall Length
WI-021-14	18.5 cm – 7 1/4

Woodson Plastic Filling Instrument Fig. 1
 Satin



Catalog #
WI-030-07

Woodson Plastic Filling Instrument Fig. 2
 Satin



Catalog #
WI-030-17



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Woodson Plastic Filling Instrument Fig. 3
Satin

Catalog #

WI-020-14



Yardley Gingival Cord Packer Fig. 1
Satin

Catalog #

WI-021-13



Zoellner Application Instrument
Satin

Catalog #

WI-021-14

Overall Length

13.5 – 5 1/2



Scalpel Handles
 Skalpgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi


Guyon Urethral Sound
 Mirror

Catalog #	Overall Length
WI-020-14	26.5 cm – 10 1/2



Sims Uterine Sound
 Satin

Catalog #	Overall Length
WI-021-13	33 cm – 13



Valleix Uterine Sound
 Satin

Catalog #	Overall Length
WI-021-14	29.5 cm – 11 1/2




Van Buren Urethral Sound
 Mirror

Catalog #	Overall Length
WI-030-07	27 cm – 10 3/4





Scalpel Handles
Skalpelligriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Antrum Curette
Satin

Catalog #	Overall Length
WI-020-14	19 cm – 7 1/2

Aspiration Uterine Curette and Coupling Handle Piece
Satin

Catalog #	Overall Length
WI-021-13	29 cm – 11 1/2



Barnhart Curette Fig. 1/2
Satin

Catalog #
WI-021-14



Barnhart Curette Fig. 5/6
Satin

Catalog #
WI-030-07



Barnhill Adenoid Curette
Satin

Catalog #	Overall Length
WI-030-17	22 cm – 8 3/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manici per bisturi

Barth Ear Scoop
 Satin



Catalog #	Overall Length
WI-020-14	18 cm – 7

Bates Curette Fig. 13 S/14S
 Satin



Catalog #
WI-021-13

Beckmann Adenoid Curette
 Satin



Catalog #
WI-021-14

Billeau Ear Loop
 Satin



Catalog #	Overall Length
WI-030-07	16 cm – 6 1/4

Blake Uterine Curette
 Satin



Catalog #	Overall Length
WI-030-17	28 cm – 11



Scalpel Handles
Skalpellgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Bruns Bone Curette
Satin

Catalog # _____
WI-020-14



Buck Ear Curette
Satin

Catalog # _____ Overall Length _____
WI-021-13 17 cm – 6 3/4



Bunting Curette Fig. 1/2
Satin

Catalog # _____
WI-021-14



Bunting Curette Fig. 5/6
Satin

Catalog # _____
WI-030-07




Caspar Bone Curette
Satin

Catalog # _____ Overall Length _____
WI-030-17 25 cm – 9 3/4



Scalpel Handles
Skalpelligriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Coakley Antrum Curette
Satin



Catalog # Overall Length
WI-020-14 17 cm – 6 3/4

Cobb Bone Curette
Satin



Catalog # Overall Length
WI-021-13 28 cm – 11

Columbia Curette Fig. 13/14
Satin



Catalog #
WI-021-14

Columbia Curette Fig. 2R/2L
Satin



Catalog #
WI-030-07

Columbia Curette Fig. 4R/4L
Satin



Catalog #
WI-030-17



Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Cureta Nasal de Halle
Satinado, Satin

Catalog #	Overall Length
WI-020-14	21 cm – 8 1/4 , 21 cm – 8 1/4



Curette Fig. 411/412
Satin

Catalog #	WI-021-13
-----------	-----------



Curette Fig. 413/414
Satin

Catalog #	WI-021-14
-----------	-----------




Curette Fig. ME-1
Satin

Catalog #	WI-030-07
-----------	-----------



Daubenspeck Bone Curette
Satin

Catalog #	Overall Length
WI-030-17	20 cm – 8



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturías
 Manches pour bistouris
 Manici per bisturi

DeMarco Furcation Curette Fig. 1
 Satin

Catalog #
 WI-020-14



DeMarco Furcation Curette Fig. 2
 Satin

Catalog #
 WI-021-13



Double Ended Bone Curette
 Satin

Catalog # Overall Length
 WI-021-14 21 cm – 8 1/4




Double Ended Molt Bone Curette
 Satin

Catalog # Overall Length
 WI-030-07 19 cm – 7 1/2



Double Ended Volkmann Bone Curette
 Satin

Catalog #
 WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Fahlbusch Micro Curette
Satin

Catalog #	Overall Length
WI-020-14	26 cm – 10 1/4



Fahlbusch Micro Scoop
Satin

Catalog #	Overall Length
WI-021-13	26 cm – 10 1/4



Faulkner Antrum Curette
Satin

Catalog #	Overall Length
WI-021-14	21 cm – 8 1/4




Ferrozell Handle Volkmann Bone Curette
Satin

Catalog #	Overall Length
WI-030-07	



Fox Lupus Curette
Satin

Catalog #	Overall Length
WI-030-17	14 cm – 5 1/2



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris
 Manici per bisturi

Goldman Curette Fig. 31/32
 Satin

Catalog #
 WI-020-14



Goldman-Fox Curette Fig. GF-2
 Satin

Catalog #
 WI-021-13



Goldman-Fox Curette Fig. GF-3
 Satin

Catalog #
 WI-021-14



Goldman-Fox Curette Fig. GF-4
 Satin

Catalog #
 WI-030-07



Gracey Curette
 Satin

Catalog #
 WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi



Greene Uterine Curette
Satin

Catalog #	Overall Length
WI-020-14	28 cm – 11



Gross Ear Scoop
Satin

Catalog #	Overall Length
WI-021-13	13 cm – 5



Hatfield Bone Curette
Satin

Catalog #	Overall Length
WI-021-14	23 cm – 9




Hebra Chalazion Curette
Satin

Catalog #	Overall Length
WI-030-07	12.5 cm – 5



Hemingway Bone Curette
Satin

Catalog #	Overall Length
WI-030-17	17 cm – 6 3/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Hoexter Periodontal Curette Fig. 1
 Satin

Catalog # Overall Length
 WI-020-14 17 cm – 6 3/4




Hoexter Periodontal Curette Fig. 2
 Satin

Catalog # Overall Length
 WI-021-13 17 cm – 6 3/4




Hoexter Periodontal Curette Fig. 3
 Satin

Catalog # Overall Length
 WI-021-14 17 cm – 6 3/4



House Micro Ear Curette
 Satin

Catalog # Overall Length
 WI-030-07 16 cm – 6 1/4



Katsch Comedone Extractor
 Satin

Catalog # Overall Length
 WI-030-17 16 cm – 6 1/4





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi


Moncorps Milium Knife
Satin

Catalog #	Overall Length
WI-020-14	14 cm – 5 1/2




Saalfeld Comedone Extractor
Satin

Catalog #	Overall Length
WI-021-13	



Unna-Vidal Comedone Extractor
Satin

Catalog #	Overall Length
WI-021-14	11.5 cm – 4 1/2



Vidal Lupus Curette
Satin

Catalog #	Overall Length
WI-030-07	13 cm – 5



Wolff Lupus Curette
Satin

Catalog #	Overall Length
WI-030-17	14 cm – 5 1/2



Scalpel Handles
Skalpgriffe
Mangos de bisturíes
bistouris
turi

Adson Periosteal Elevator
Satin

Catalog # _____
WI-020-14



Alderdyce Periosteal Elevator
Satin

Catalog # _____ Overall Length _____
WI-021-13 21 cm – 8 1/4



Alexander Periosteal Elevator
Satin

Catalog # _____ Overall Length _____
WI-021-14 21 cm – 8 1/4



Allen Periosteal Raspatory / Elevator
Satin

Catalog # _____
WI-030-07



Allen Periosteal Raspatory / Elevator Fig. 9A
Satin

Catalog # _____
WI-030-17





Scalpel Handles
Skalpelligriffe
Mangos de bisturías
Manches pour bistouris
Manici per bisturi



Apexo Root Elevator Fig. 301
Satin

Catalog # _____
WI-020-14



Apexo Root Elevator Fig. 302
Satin

Catalog # _____
WI-021-13



Apexo Root Elevator Fig. 303
Satin

Catalog # _____ Overall Length
WI-021-14 11.5 cm – 4 1/2



Apical Root Tip Pick
Satin

Catalog # _____
WI-030-07



Austin Attic Dissector
Satin

Catalog # _____ Overall Length
WI-030-17 15.5 cm – 6



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris
 Manici per bisturi

Ballanger-Hajek Septum Elevator
 Satin

Catalog #	Overall Length
WI-020-14	19 cm – 7 1/2



Double Ended, 6.0 mm, 6.0 mm

Beales Elevator
 Satin

Catalog #	Overall Length
WI-021-13	16.5 cm – 6 1/2



Large
 Small


Bennet Periosteal Raspatory / Elevator Fig. 12 S
 Satin

Catalog #
WI-021-14



Bennett Bone Lever (Retractor)
 Satin

Catalog #
WI-030-07



60 mm
 43 mm

Benque Periosteal Raspatory / Elevator
 Satin

Catalog #
WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Bernard Root Elevator (Syndesmotome)
Satin

Catalog #

WI-020-14



Binder Flap Elevator
Satin

Catalog #

WI-021-13

Overall Length

13 cm – 5



Blumenthal Conjunctiva Dissector
Satin

Catalog #

WI-021-14

Overall Length

13 cm – 5



Boies Nasal Fracture Elevator
Satin

Catalog #

WI-030-07

Overall Length

19.5 cm – 7 3/4




Bone Lever (Retractor)
Satin

Catalog #

WI-030-17

Overall Length


22 cm – 8 3/4



Scalpel Handles
 Skalpelligriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Bradley Root Elevator Fig. 44C
 Satin

Catalog #
 WI-020-14



Bristow Bone Lever (Retractor)
 Satin

Catalog # Overall Length
 WI-021-13 23 cm – 9



Buchs Periosteal Elevator
 Satin

Catalog # Overall Length
 WI-021-14 19 cm – 7 1/2



Buratto Double Ended Lasik Flap Protector
 Satin

Catalog # Overall Length
 WI-030-07 11.5 cm – 4 1/2



Buratto Lasik Flap Protector
 Satin

Catalog # Overall Length
 WI-030-17 12.5 cm – 5





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Buser Periosteal Raspatory / Elevator
Satin

Catalog #

WI-020-14



Carpenter Tonsil Dissector
Satin

Catalog #

WI-021-13



Caspar Micro Dissector
Satin

Catalog #

WI-021-14

Overall Length

24 cm – 9 1/2



Caspar Periosteal Elevator
Satin

Catalog #

WI-030-07



Chandler Bone Lever (Retractor)
Satin

Catalog #

WI-030-17



Scalpel Handles
 Skalpgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Chompret Root Elevator (Syndesmotome)
 Satin



Catalog # _____
 WI-020-14

Cobb Periosteal Elevator
 Satin




Catalog #	Overall Length
WI-021-13	28 cm – 11

Cogswell Root Elevator Fig. A
 Satin



Catalog # _____
 WI-021-14

Corneal Dissector
 Satin



Catalog #	Overall Length
WI-030-07	12 cm – 4 3/4

Cottle Periosteal Elevator
 Satin



Catalog # _____
 WI-030-17



Scalpel Handles
Skalpelligriffe
Mangos de bisturís
Manches pour bistouris
Manici per bisturi



Cottle Periosteal Elevator
Satin

Catalog # Overall Length
WI-020-14 19.5 cm – 7 3/4



Cottle Septum Elevator
Satin

Catalog # Overall Length
WI-021-13 23.5 cm – 9 1/4



Coupland Root Elevator
Satin

Catalog #
WI-021-14



Crane Pick Root Elevator Fig. 8
Satin

Catalog #
WI-030-07



Crego Bone Lever (Elevator)
Satin

Catalog # Overall Length
WI-030-17 19 cm – 7 1/2



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Cryer Root Elevator Fig. 21
Satin

Catalog #

WI-020-14



Cryer Root Elevator Fig. 22
Satin

Catalog #

WI-021-13



Cryer Root Elevator Fig. 25
Satin

Catalog #

WI-021-14



Cryer Root Elevator Fig. 26
Satin

Catalog #

WI-030-07



Cryer Root Elevator Fig. 31
Satin

Catalog #

WI-030-17





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Cryer Root Elevator Fig. 32
Satin

Catalog # _____
WI-020-14



Cryer Root Elevator Fig. 45
Satin

Catalog # _____
WI-021-13



Cryer Root Elevator Fig. 46
Satin

Catalog # _____
WI-021-14



Curtis Root Elevator Fig. 51
Satin

Catalog # _____
WI-030-07



Curtis Root Elevator Fig. 52
Satin

Catalog # _____
WI-030-17



Scalpel Handles
 Skalpgriffe
 Mangos de bisturíes
 pour bistouris
 or bisturi

Cushing Periosteal Elevator
 Satin

Catalog # _____
 WI-020-14



Dastoor Lacrimal Sac Dissector
 Satin

Catalog # _____ Overall Length
 WI-021-13 13 cm – 5



Fickling Nail Elevator
 Satin

Catalog # _____
 WI-021-14



Mackay Nail Elevator
 Satin

Catalog # _____ Overall Length
 WI-030-07 16.5 cm – 6 1/2



Nail Elevator
 Satin

Catalog # _____ Overall Length
 WI-030-17 16 cm – 6 1/4





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Aufricht Glabella Rasp
Satin

Catalog # Overall Length
WI-020-14 21 cm – 8 1/4



Baka Nail File
Satin

Catalog # Overall Length
WI-021-13 12 cm – 4 3/4



Baka Nail File
Satin

Catalog # Overall Length
WI-021-14 12 cm – 4 3/4




Bier Amputation Saw
Satin

Catalog # Overall Length
WI-030-07 40.5 cm – 16



Blade for Charriere Amputation Saw
Satin


Catalog #
WI-030-17



Scalpel Handles
Skalpelligriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Bone File
Satin

Catalog # Overall Length
WI-020-14 19 cm – 7 1/2



Bone Rasp
Satin

Catalog #
WI-021-13



Buck Periodontal Files Fig. 11/12
Satin

Catalog #
WI-021-14



Charriere Amputation Saw
Mirror

Catalog # Overall Length
WI-030-07 35 cm – 13 3/4



With Saw Blade

Charriere Metacarpal Saw
Mirror

Catalog # Overall Length
WI-030-17 26.5 cm – 10 1/2



With Saw Blade



Scalpel Handles
Skalpelligriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Cottle Nasal Rasp
Satin




Catalog #	Overall Length
WI-020-14	20.5 cm – 8

Diamond Deb Foot Dresser
Mirror



Catalog #	Overall Length
WI-021-13	20 cm – 8

Dingman Otoprader
Satin



Catalog #	Overall Length
WI-021-14	15 cm – 6

Engel Plaster Saw
Satin



Catalog #	Overall Length
WI-030-07	15.5 cm – 6

Finger Ring Saw
Satin



Catalog #	Overall Length
WI-030-17	



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manici per bisturi


Finger Ring Saw Blade
 Satin

Catalog # _____
 WI-020-14




Flat Bone File
 Satin

Catalog # _____
 WI-021-13




Fomon Bone File
 Satin

Catalog # _____ Overall Length _____
 WI-021-14 21 cm – 8 1/4



Fomon Nasal Rasp
 Mirror

Catalog # _____ Overall Length _____
 WI-030-07 20.5 cm – 8



Foot Dresser
 Mirror

Catalog # _____ Overall Length _____
 WI-030-17 20 cm – 8





Scalpel Handles
Skalpellgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Gallaher Antrum Rasp
Satin

Catalog #	Overall Length
WI-020-14	17 cm – 6 3/4



Gigli Saw Handle
Mirror

Catalog #
WI-021-13



Gigli Wire Saw
Satin

Catalog #
WI-021-14




Green Emry Dressor
Satin

Catalog #	Overall Length
WI-030-07	18 cm – 7



Joseph Bone File
Satin

Catalog #	Overall Length
WI-030-17	16 cm – 6 1/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris
 Manici per bisturi

Langenbeck Amputation Saw
 Satin




Catalog #	Overall Length
WI-020-14	24.5 cm – 9 3/4

Lewis Nasal Rasp
 Satin



Catalog #	Overall Length
WI-021-13	19 cm – 7 1/2

Maltz Nasal Rasp
 Satin



Catalog #	Overall Length
WI-021-14	17 cm – 6 3/4

Martin Finger Ring Saw
 Satin



Catalog #	Overall Length
WI-030-07	18.5 cm – 7 1/4

Medium Tip Nail File
 Satin




Catalog #
WI-030-17



Scalpel Handles
Skalpgriffe
Mangos de bisturís
Manches pour bistouris
Manici per bisturi


Miller Bone File
Satin

Catalog #	Overall Length
WI-020-14	19 cm – 7 1/2



Miller-Colburn Bone File
Satin

Catalog #	Overall Length
WI-021-13	18 cm – 7



Nail File
Satin

Catalog #
WI-021-14



Medium Tip, Fine / Cross Cut

Parkes Bone File
Satin

Catalog #	Overall Length
WI-030-07	19.5 cm – 7 3/4



Putti Bone Rasp
Satin

Catalog #
WI-030-17



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris
 Manici per bisturi

Rhein Amalgam File Fig. 31/32
 Satin



Catalog # _____
 WI-020-14

Satterlee Saw Blade
 Mirror



Catalog # Overall Length
 WI-021-13 29 cm – 11 1/2

Schluger Periodontal Files Fig. 9/10
 Satin



Catalog # _____
 WI-021-14

Semi Round Bone File
 Satin



Catalog # _____
 WI-030-07

Sugarman Periodontal Files Fig. 3S/4S
 Satin



Catalog # _____
 WI-030-17



Scalpel Handles
Skalpelligriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Walton Finger Ring Saw
Satin

Catalog #	Overall Length
WI-020-14	16 cm – 6 1/4




Waqar Intl OsciPower™ Oscillating Plaster Saw

Catalog #	
WI-021-13	



Webster Glabella Nasal Rasp
Satin

Catalog #	Overall Length
WI-021-14	21 cm – 8 1/4




Wedelstaedt Amalgam File Fig. 1/4
Satin

Catalog #	
WI-030-07	



Wiener Nasal Rasp
Satin

Catalog #	Overall Length
WI-030-17	18 cm – 7



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Acland Clip Applying Forcep
 Satin




Catalog #	Overall Length
WI-020-14	14.5 cm – 5 3/4

Adair Tenaculum Forcep
 Satin




Catalog #	Overall Length
WI-021-13	19 cm – 7 1/2

Adlercreutz Dissecting Forcep
 Satin



Catalog #
WI-021-14

Adson Dissecting Forcep
 Satin



Catalog #
WI-030-07

Adson Ear Forceps
 Satin




Catalog #	Overall Length
WI-030-17	17.5 cm – 7



Scalpel Handles
Skalpgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Adson-Brown Dissecting Forcep
Satin

Catalog #	Overall Length
WI-020-14	12 cm – 5 , 12.5 cm – 4 3/4 , 15 cm – 6



Aiach Columella Suturing Forceps
Satin

Catalog #	Overall Length
WI-021-13	15.5 cm – 6



Akahoshi Phaco PreChopper
Satin

Catalog #	Overall Length
WI-021-14	11 cm – 4 1/4



Alfonso Nucleus Grasping Forceps
Satin

Catalog #	Overall Length
WI-030-07	10.5 cm – 4 1/4



American Pattern Tooth Extraction Forceps (Child) Fig. 150S
Mirror

Catalog #
WI-030-17



Scalpel Handles
Skalpellgriffe
Mangos de bisturís
Manches pour bistouris

American Pattern Tooth Extraction Forceps (Child) Fig. 151S
Mirror

Catalog # _____

WI-020-14



American Pattern Tooth Extraction Forceps Fig. 1
Mirror

Catalog # _____

WI-021-13



American Pattern Tooth Extraction Forceps Fig. 10S
Mirror

Catalog # _____

WI-021-14



American Pattern Tooth Extraction Forceps Fig. 150AS
Mirror

Catalog # _____

WI-030-07



American Pattern Tooth Extraction Forceps Fig. 151AS
Mirror

Catalog # _____

WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



American Pattern Tooth Extraction Forceps Fig. 16
Mirror

Catalog # _____

WI-020-14



American Pattern Tooth Extraction Forceps Fig. 16S
Mirror

Catalog # _____

WI-021-13



American Pattern Tooth Extraction Forceps Fig. 17
Mirror

Catalog # _____

WI-021-14



American Pattern Tooth Extraction Forceps Fig. 1S
Mirror

Catalog # _____

WI-030-07



American Pattern Tooth Extraction Forceps Fig. 203
Mirror

Catalog # _____

WI-030-17



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manici per bisturi

American Pattern Tooth Extraction Forceps Fig. 210H
 Mirror

Catalog # _____
 WI-020-14



American Pattern Tooth Extraction Forceps Fig. 210S
 Mirror

Catalog # _____
 WI-021-13



American Pattern Tooth Extraction Forceps Fig. 217
 Mirror

Catalog # _____
 WI-021-14



American Pattern Tooth Extraction Forceps Fig. 222
 Mirror


Catalog # _____
 WI-030-07



American Pattern Tooth Extraction Forceps Fig. 286
 Mirror

Catalog # _____
 WI-030-17





Scalpel Handles
Skalpelli-griffe
Mangos de bisturís
Manches pour bistouris
Manici per bisturi

American Pattern Tooth Extraction Forceps Fig. 39
Mirror

Catalog # _____
WI-020-14



American Pattern Tooth Extraction Forceps Fig. 40
Mirror

Catalog # _____
WI-021-13



American Pattern Tooth Extraction Forceps Fig. 53L
Mirror

Catalog # _____
WI-021-14



American Pattern Tooth Extraction Forceps Fig. 53R
Mirror

Catalog # _____
WI-030-07



American Pattern Tooth Extraction Forceps Fig. 62
Mirror

Catalog # _____
WI-030-17



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

American Pattern Tooth Extraction Forceps Fig. 65
 Mirror

Catalog #
 WI-020-14



Anvil Nail Extracting Forceps
 Satin

Catalog # Overall Length
 WI-021-13 12.5 cm – 5



Arruga Capsule Forcep
 Satin

Catalog # Overall Length
 WI-021-14 11 cm – 4 1/4




Asch Septum Straightning Forceps
 Satin

Catalog # Overall Length
 WI-030-07 22.5 cm – 8 3/4



Austin Micro Suturing Forcep
 Satin

Catalog # Overall Length
 WI-030-17 14 cm – 5 1/2





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Baird Chalazion Forcep
Satin

Catalog #	Overall Length
WI-020-14	9.5 cm – 3 3/4



Bechert Lens Holding Forcep
Satin

Catalog #	Overall Length
WI-021-13	12.5 cm – 5



Beer Cilia Forcep
Satin

Catalog #	Overall Length
WI-021-14	9 cm – 3 1/2




Bergh Cilia Forcep
Satin

Catalog #	Overall Length
WI-030-07	8.5 cm – 3 1/4



Berke Ptosis Forcep
Satin

Catalog #	Overall Length
WI-030-17	11 cm – 4 1/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manîni per bisturi

Bill's Traction Handle
 Satin

Catalog # _____
 WI-020-14



Bircher-Ganske Cartilage Seizing Forcep
 Satin

Catalog # _____
 WI-021-13



Bishop-Harmon Dressing Forcep
 Satin

Catalog #	Overall Length
WI-021-14	8.5 cm – 3 1/4



Cross Serrated Jaws

Bohler Bone Cutting Forcep
 Satin

Catalog #	Overall Length
WI-030-07	14.5 cm – 5 3/4



Compound Action

Bonaccollo Conjunctival Forcep
 Satin

Catalog #	Overall Length
WI-030-17	10.5 cm – 4





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi


Bone Holding & Repositioning Forcep
Satin

Catalog #	Overall Length
WI-020-14	14 cm – 5 1/2



Bone Holding & Repositioning Forcep with Thread Fixation
Satin

Catalog #
WI-021-13



Bonn Corneal Forcep
Satin

Catalog #
WI-021-14



Brand Tendon Tunnelling Forcep
Satin

Catalog #
WI-030-07




Braun Cranioclast
Satin

Catalog #	Overall Length
WI-030-17	42 cm – 16 1/2




Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Brophy Dissecting Forcep
 Satin




Catalog #	Overall Length
WI-020-14	20 cm – 8

Brophy Dressing Forcep
 Satin



Catalog #	Overall Length
WI-021-13	20 cm – 8

Buck Foreign Body Forceps
 Satin



Catalog #	Shaft Length
WI-021-14	11.5 cm – 4 1/2

Chevalier-Jackson Grasping Forceps
 Satin



Catalog #
WI-030-07

Fechtner Conjunctival Forcep
 Satin



Very Delicate Ring Stopped Tips With Tying Platform

Catalog #	Overall Length
WI-030-17	10.5 cm – 4



Scalpel Handles
Skalpelligriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Lerche Dissecting Forceps
Satin

Catalog #	Overall Length
WI-020-14	15 cm – 6



Mod. USA Dissecting Forcep
Satin

Catalog #	Overall Length
WI-021-13	



Radolf Nail Extracting Forceps
Satin

Catalog #	Overall Length
WI-021-14	14 cm – 5 1/2




Slender Pattern Dissecting Forcep
Satin

Catalog #	Overall Length
WI-030-07	




Slender Pattern Dressing Forcep
Satin

Catalog #	Overall Length
WI-030-17	



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Standard Dissecting Forcep
 Satin



Catalog # _____

WI-020-14


Standard Dressing Forcep
 Satin



Catalog # _____

WI-021-13

Stille Dissecting Forcep
 Satin



Catalog # _____

Overall Length _____

WI-021-14

15 cm – 6

Stone Dissecting Forcep
 Satin



Catalog # _____

Overall Length _____

WI-030-07

15 cm – 6

Tendon Interlacing Forceps
 Satin



Catalog # _____


WI-030-17



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Tunneling Forcep
Satin

Catalog #	Overall Length
WI-020-14	21 cm – 8 1/4



Acland Double Micro Vessel Clamp
Satin

Catalog #
WI-021-13



Acland Micro Vessel Clamp
Satin

Catalog #
WI-021-14




Adair Tenaculum Forcep
Satin

Catalog #	Overall Length
WI-030-07	19 cm – 7 1/2



Adson-Baby Dissecting and Ligature Forcep
Satin

Catalog #
WI-030-17



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Allen Intestinal Clamp
 Satin

Catalog #	Overall Length
WI-020-14	16 cm – 6 1/4



Allis Atrauma Intestinal and Tissue Grasping Forcep
 Satin

Catalog #
WI-021-13



Allis Intestinal and Tissue Grasping Forcep
 Satin

Catalog #
WI-021-14



Allis Intestinal and Tissue Grasping Forcep
 Satin

Catalog #	Overall Length
WI-030-07	20 cm – 8




Allis-Adair Intestinal and Tissue Grasping Forcep
 Satin

Catalog #	Overall Length
WI-030-17	16 cm – 6 1/4





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Allis-Baby Intestinal and Tissue Grasping Forcep
Satin

Catalog #	Overall Length
WI-020-14	13 cm – 5



Atrauma Organ Grasping / Haemorrhoidal Forcep
Satin

Catalog #	Overall Length
WI-021-13	



Babcock Atrauma Intestinal and Tissue Grasping Forcep
Satin

Catalog #	Overall Length
WI-021-14	




Babcock Intestinal and Tissue Grasping Forcep
Satin

Catalog #	Overall Length
WI-030-07	



Baby-Satinsky Anastomosis Clamps
Satin

Catalog #	Overall Length
WI-030-17	15 cm – 6



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
lanches pour bistouris
lanici per bisturi

Backhaus Towel Clamp
Satin

Catalog # _____
WI-020-14



Bainbridge Artery Forcep
Satin

Catalog # _____
WI-021-13




Bainbridge Atrauma Multipurpose Vascular Clamps
Satin

Catalog # _____
WI-021-14




Barre Dissecting and Ligature Forcep
Satin

Catalog #	Overall Length
WI-030-07	28 cm – 11



Bengolea Artery Forcep
Satin

Catalog # _____
WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi




Bergmann Sponge Holding Forcep
Satin

Catalog #	Overall Length
WI-020-14	24 cm – 9 1/2




Berkeley-Bonney Hysterectomy Forcep
Satin

Catalog #	Overall Length
WI-021-13	20 cm – 8



Bernhard Towel Clamp
Satin

Catalog #	Overall Length
WI-021-14	16.5 cm – 6 1/2




Biemer Clip Applying Forcep
Satin

Catalog #	Overall Length
WI-030-07	14.5 cm – 5 3/4



Biemer Vessel Clip
Satin

Catalog #	Overall Length
WI-030-17	




Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manici per bisturi



Biemer-Muller Approximator
 Satin

Catalog # _____
 WI-020-14

Birkett Artery Forcep
 Satin




Catalog # Overall Length
 WI-021-13 18.5 cm – 7 1/4

Birkett Artery Forcep
 Satin



Catalog # Overall Length
 WI-021-14 18.5 cm – 7 1/4

Blake Gall Stone Forcep
 Satin



Catalog # Overall Length
 WI-030-07 21 cm – 8 1/4

Blohmke Tonsil Grasping Forceps
 Satin



Catalog # Overall Length
 WI-030-17 20 cm – 8



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Bonney Sponge Holding Forcep
Satin

Catalog #	Overall Length
WI-020-14	24 cm – 9 1/2



Boys-Allis Intestinal and Tissue Grasping Forcep
Satin

Catalog #	Overall Length
WI-021-13	15.5 cm – 6



Bozeman Sponge Holding Forcep
Satin

Catalog #	Overall Length
WI-021-14	



Bozemann Knot Tier
Satin

Catalog #	Overall Length
WI-030-07	21 cm – 8 1/4



Braun Tenaculum Forcep
Satin


Catalog #	Overall Length
WI-030-17	



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris
 Manici per bisturi

Bridge Artery Forcep
 Satin

Catalog #	Overall Length
WI-020-14	27.5 cm – 10 3/4



Brompton Tube Milking Rolle
 Blue Anodized

Catalog #	Overall Length
WI-021-13	



Brunner Atrauma Sigmoid Anastomosis Clamp
 Satin

Catalog #	Overall Length
WI-021-14	27 cm – 10 3/4



Brunner Ligature Conductor
 Satin


Catalog #	Overall Length
WI-030-07	18 cm – 7



Ferrozell Solid Handle

Brunner Sigmoid Anastomosis Clamp
 Satin

Catalog #	Overall Length
WI-030-17	24 cm – 9 1/2





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi




Bryant Dressing Forcep
Satin

Catalog #	Overall Length
WI-020-14	13 cm – 5




Burke Hysterectomy Forcep

Catalog #	Overall Length
WI-021-13	25 cm – 9 3/4



Cairns Artery Forcep
Satin

Catalog #	Overall Length
WI-021-14	



Carrel Haemostatic Forcep
Satin

Catalog #	Overall Length
WI-030-07	8.5 cm – 3 1/4



Caspar Ligature Adjuster
Satin

Catalog #	Overall Length
WI-030-17	23 cm – 9



Bayonet Shaped



Scalpel Handles
Skalpgriffe
Mangos de bisturís
Manches pour bistouris
Manici per bisturi

Chaput Tissue Grasping Forcep
Satin

Catalog #	Overall Length
WI-020-14	12.5 cm – 5



Cheron Dressing Forcep
Satin

Catalog #	Overall Length
WI-021-13	23.5 cm – 9 1/4




Childe Approximation Forcep
Satin

Catalog #	Overall Length
WI-021-14	17.5 cm – 7



Tunneling Forcep
Satin

Catalog #	Overall Length
WI-030-07	21 cm – 8 1/4



Denture Cup
Satin

Catalog #	Overall Length
WI-030-17	Ø 75 x 50 mm



With Lid



Scalpel Handles
Skalpellgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Dreassing Jar
Satin

Catalog # _____
WI-020-14



Female Urinal
Satin

Catalog # _____
WI-021-13



Gallipot / Mixing Bowl
Satin

Catalog # _____
WI-021-14



HW-021-14



HW-021-11

Graduated Measure
Satin

Catalog # _____
WI-030-07



HW-070-01



HW-072-01

Graduated Triangular Basin
Satin

Catalog # _____ Overall Length
WI-030-17 18 cm – 7



With Handle



Without Handle



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Instrument Box
Satin

Catalog # _____
WI-020-14



Instrument Box with Wire Lock
Satin

Catalog # _____ Overall Length _____
WI-021-13 205 x 105 x 35 mm



Instrument Jar
Satin

Catalog # _____
WI-021-14



Instrument Tray
Satin

Catalog # _____
WI-030-07



Kidney Dish
Satin

Catalog # _____
WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi



Male Urinal
Satin

Catalog # _____
WI-020-14



Medicine Cup
Satin

Catalog # _____
WI-021-13



Medicine Spoon
Satin

Catalog # Overall Length
WI-021-14 20 ccm



Metal Cup
Satin

Catalog # Overall Length
WI-030-07 Ø 90 x 90 mm



With Lid

Methot Cotton Dispenser
Satin

Catalog # Overall Length
WI-030-17 Ø 50 x 75 mm



Without Lid



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistour
 Manici per bisturi



HI-020-04



HI-020-11

Mixing Bowl
 Satin

Catalog # _____
 WI-020-14

Perfection Type Bed Pan
 Satin

Catalog # _____
 WI-021-13



Petri Dish
 Satin

Catalog #	Overall Length
WI-021-14	Ø 100 x 20 mm



Placenta Pan with Cover
 Satin

Catalog # _____
 WI-030-07



Round Type Bed Pan
 Satin

Catalog # _____
 WI-030-17





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Spirit Lamp
Satin

Catalog #	Overall Length
WI-020-14	120 ccm



Sputum Cup
Satin

Catalog #	Overall Length
WI-021-13	Ø 75 x 80 mm



Stainless Steel Basin
Satin

Catalog #	Overall Length
WI-021-14	



Tablet Tray
Satin

Catalog #	Overall Length
WI-030-07	Ø 75 x 12 mm



Universal Instrument Tray
Satin

Catalog #	Overall Length
WI-030-17	



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Alexander Bone Chisel
 Dual Satin & Mirror

Catalog #	Overall Length
WI-020-14	17.5 cm – 7



Alexander Bone Gouge
 Dual Satin & Mirror

Catalog #	Overall Length
WI-021-13	17.5 cm – 7



Bauer Osteotome
 Dual Satin & Mirror

Catalog #	Overall Length
WI-021-14	20.5 cm – 8



Bone Chisel
 Dual Satin & Mirror

Catalog #	Overall Length
WI-030-07	13.5 cm – 5 1/4



Bone Gouge
 Dual Satin & Mirror

Catalog #	Overall Length
WI-030-17	13.5 cm – 5 1/4





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Bone Gouge for Lamina
Dual Satin & Mirror

Catalog # Overall Length
WI-020-14 23 cm – 9 1/4



Bone Mallet
Satin

Catalog #
WI-021-13



Ø 30.0 mm, Ø 42.0 mm, 285 Grams

Bone Osteotome
Dual Satin & Mirror

Catalog # Overall Length
WI-021-14 13.5 cm – 5 1/4






Bone Spreading Osteotome
Satin

Catalog #
WI-030-07



Bone Tamper
Satin

Catalog # Overall Length
WI-030-17 15.5 cm – 6



Scalpel Handles
Skalpellgriffe
Mangos de bisturios
Manches pour bistouris

Buck Percussion Reflex Hammer
Satin

Catalog # Overall Length
WI-020-14 18.5 cm – 7 1/4




Bunnell Hand Drill
Satin

Catalog #
WI-021-13



Cinzel de Cottle
Dual Satin & Mirror, Dual Satin & Mirror

Catalog # Overall Length
WI-021-14 18.5 cm – 7 1/4 , 18.5 cm – 7 1/4



Cinzel Estándar
Dual Satin & Mirror, Dual Satin & Mirror

Catalog # Overall Length
WI-030-07 13.5 cm – 5 1/4 , 13.5 cm – 5 1/4




Cinelli Bone Osteotome
Dual Satin & Mirror

Catalog # Overall Length
WI-030-17 16 cm – 6 1/4





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Cloward Bone Gouge
Dual Satin & Mirror

Catalog #	Overall Length
WI-020-14	24.5 cm – 9 3/4



Cloward Bone Osteotome
Dual Satin & Mirror

Catalog #	Overall Length
WI-021-13	24.5 cm – 9 3/4



Collin Bone Mallet
Satin


Catalog #	Overall Length
WI-021-14	20 cm – 8



Ø 30.0 mm, 220 Grams


Cottle Bone Chisel
Dual Satin & Mirror

Catalog #	Overall Length
WI-030-07	18.5 cm – 7 1/4



Cottle Bone Osteotome
Dual Satin & Mirror

Catalog #	Overall Length
WI-030-17	18.5 cm – 7 1/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Dahmen Bone Osteotome
 Dual Satin & Mirror

Catalog #	Overall Length
WI-020-14	30 cm – 12



Ferrozell Handle Bone Chisel
 Dual Satin & Mirror

Catalog #	Overall Length
WI-021-13	



Ferrozell Handle Bone Tamper
 Satin

Catalog #	Overall Length
WI-021-14	20 cm – 8



Gardner Bone Chisel
 Dual Satin & Mirror

Catalog #	Overall Length
WI-030-07	17 cm – 6 3/4



Hibbs Bone Chisel
 Dual Satin & Mirror

Catalog #	Overall Length
WI-030-17	24.5 cm – 9 3/4





Scalpel Handles
Skalpellsgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Hibbs Bone Osteotome
Dual Satin & Mirror

Catalog #	Overall Length
WI-020-14	24.5 cm – 9 3/4




Hoke Bone Osteotome
Dual Satin & Mirror

Catalog #	Overall Length
WI-021-13	17 cm – 6 3/4



Lambotte Bone Osteotome
Dual Satin & Mirror

Catalog #	Overall Length
WI-021-14	24 cm – 9 1/2



Lexer Bone Chisel
Dual Satin & Mirror

Catalog #	Overall Length
WI-030-07	20 cm – 8



Lexer Bone Osteotome
Dual Satin & Mirror

Catalog #	Overall Length
WI-030-17	



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Lucas Bone Chisel
 Dual Satin & Mirror

Catalog # Overall Length
 WI-020-14 17 cm – 6 3/4




Lucas Bone Osteotome
 Dual Satin & Mirror

Catalog # Overall Length
 WI-021-13 17 cm – 6 3/4



Mini-Lambotte Bone Osteotome
 Dual Satin & Mirror

Catalog #
 WI-021-14



Neuro Tuning Fork
 Satin

Catalog #
 WI-030-07



With Hammer, Pinwheel, Brush and Needle

Partsch Bone Chisel
 Dual Satin & Mirror

Catalog #
 WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Partsch Bone Gouge
Dual Satin & Mirror

Catalog # _____
WI-020-14



Post Mortem Hammer with Hook Handle
Satin

Catalog # _____ Overall Length _____
WI-021-13 24 cm – 9 1/2



Rabiner Percussion Hammer
Satin

Catalog # _____ Overall Length _____
WI-021-14 23 cm – 9



Rachiotome
Satin

Catalog # _____ Overall Length _____
WI-030-07 22 cm – 8 3/4



Skull Breaker & Chisel
Satin


Catalog # _____ Overall Length _____
WI-030-17 14 cm – 5 1/2



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris

Smith-Peterson Bone Osteotome
 Dual Satin & Mirror

Catalog #	Overall Length
WI-020-14	20.5 cm – 8



Snyder OSU Bone Gouge
 Dual Satin & Mirror

Catalog #	Overall Length
WI-021-13	16.5 cm – 6 1/2



Stille Bone Chisel
 Dual Satin & Mirror

Catalog #	Overall Length
WI-021-14	20 cm – 8



Stille Bone Osteotome
 Dual Satin & Mirror

Catalog #	Overall Length
WI-030-07	20 cm – 8



Taylor Percussion Hammer
 Satin

Catalog #	Overall Length
WI-030-17	20 cm – 8





Scalpel Handles
Skalpellgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi




Virchow Skull Breaker
Satin

Catalog # _____
WI-020-14




Wagner Bone Gouge
Dual Satin & Mirror

Catalog # _____ Overall Length _____
WI-021-13 34 cm – 13 1/2




Walter Bone Chisel
Dual Satin & Mirror

Catalog # _____ Overall Length _____
WI-021-14 19 cm – 7 1/2



Walter Bone Osteotome
Dual Satin & Mirror

Catalog # _____ Overall Length _____
WI-030-07 19 cm – 7 1/2



Wartenberg Sensibility Instrument
Satin

Catalog # _____ Overall Length _____
WI-030-17 19 cm – 7 1/2



Scalpel Handles
 Skalpelligriffe
 Mangos de bisturías
 Manches pour bistouris
 Manici per bisturi

Anillo de Estabilización de Gimbel
 Satinado, Satin



Catalog #	Overall Length
WI-020-14	11.5 cm – 4 1/2 , 11.5 cm – 4 1/2

Anis Radial Marker
 Satin




Catalog #	Overall Length
WI-021-13	10 cm – 4

Anis Suture Placement Marker
 Satin



Catalog #	Overall Length
WI-021-14	10 cm – 4

Boley Gauge
 Satin



Catalog #	
WI-030-07	

Bores Optic Zone Marker
 Satin



Catalog #	Overall Length
WI-030-17	10.5 cm – 4 1/4



Scalpel Handles
Skalpelligriffe
Mangos de bisturias
Manches pour bistouris
Manici per bisturi



Chayet Corneal Marker
Satin

Catalog #	Overall Length
WI-020-14	11.5 cm – 4 1/2



Fechtner Trabeculectomy Marker
Satin

Catalog #	Overall Length
WI-021-13	9 cm – 3 1/2



Fine-Thornton Swivel Fixation Ring
Satin

Catalog #	Overall Length
WI-021-14	12 cm – 4 3/4



Grene Visual Axis Marker
Satin

Catalog #	Overall Length
WI-030-07	11.5 cm – 4 1/2



Helveston Scleral Ruler
Satin

Catalog #	Overall Length
WI-030-17	12 cm – 4 3/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Hofmann-Thornton Globe Fixation Ring
 Satin



Catalog #	Overall Length
WI-020-14	11 cm – 4 1/4

House Micro Measuring Rod
 Satin



Catalog #	Overall Length
WI-021-13	15.5 cm – 6

Iwanson Caliper
 Satin



Catalog #	
WI-021-14	

Jameson Caliper
 Satin



Catalog #	
WI-030-07	

Kellen Capsulorhexis Marker
 Satin



Catalog #	Overall Length
WI-030-17	11.5 cm – 4 1/2



Scalpel Handles
Skalpellgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Maloney Intra-Operative Keratometer
Satin

Catalog # _____
WI-020-14



6 mm Diameter Marking Ring
With 2.5 mm Extensions

McDonald Optic Zone Marker
Satin

Catalog # _____ Overall Length _____
WI-021-13 11 cm – 4 1/4



Mendez Degree Gauge
Satin

Catalog # _____ Overall Length _____
WI-021-14 11.5 cm – 4 1/2



Metal Ruler
Satin

Catalog # _____
WI-030-07



Graduated in mm and inches
150 mm, 200 mm, 300 mm, 500 mm

Schiötz Tonometer
Satin


Catalog # _____
WI-030-17



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Stahl Corneal Caliper
 Satin

Catalog # _____
 WI-020-14



Thornton Globe Fixation Ring
 Satin

Catalog # _____ Overall Length
 WI-021-13 10 cm – 4



Wartenberg Sensibility Instrument
 Satin

Catalog # _____ Overall Length
 WI-021-14 19 cm – 7 1/2



Abraham Tonsil Knife
 Satin

Catalog # _____
 WI-030-07



Andrew Carver Fig. 1
 Satin

Catalog # _____
 WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Andrew Carver Fig. 2
Satin

Catalog #

WI-020-14



Andrew Carver Fig. 3
Satin

Catalog #

WI-021-13



Ayre Conization Knife
Satin

Catalog #

WI-021-14

Overall Length

23.5 cm – 9 1/4



Ballanger Swivel Knife
Satin

Catalog #

WI-030-07



Bayonet Shaped
Straight

Barbara Micro Ear Needle
Satin

Catalog #

WI-030-17

Overall Length

16 cm – 6 1/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris

Barbara Micro Ear Pick
 Satin

Catalog #	Overall Length
WI-020-14	16 cm – 6 1/4



Bard-Parker Scalpel Handle No. 3
 Satin

Catalog #	Overall Length
WI-021-13	12 cm – 4 3/4



Barkan Goniotomy Knife
 Satin

Catalog #	Overall Length
WI-021-14	11.5 cm – 4 1/2



Barraquer Blade Breaker
 Satin

Catalog #	Overall Length
WI-030-07	11 cm – 4 1/4



Beaver Scalpel Handle
 Satin

Catalog #	Overall Length
WI-030-17	



For Micro Scalpel Blades



Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Bisturí para Colgajo de Plester
Satinado, Satin



Catalog #	Overall Length
WI-020-14	15 cm – 6 , 15 cm – 6

Black Carving Knife Fig. 7
Satin



Catalog #
WI-021-13

Buck Periodontal Knife Fig. 3/4
Satin



Catalog #
WI-021-14

Buck Periodontal Knife Fig. 5/6
Satin



Catalog #
WI-030-07

Carpenter Tonsil Knife
Satin



Catalog #
WI-030-17



Scalpel Handles
Skalpelli-griffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Carver Fig. 104
Satin



Catalog #

WI-020-14

Carving Knife Fig. 14L
Satin



Catalog #

WI-021-13

Carving Knife Fig. 29
Satin



Catalog #

WI-021-14

Carving Knife Fig. 2S
Satin



Catalog #

WI-030-07

Carving Knife Fig. 6
Satin



Catalog #

WI-030-17



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Caspar Scalpel Handle
Satin

Catalog #	Overall Length
WI-020-14	18 m – 7



Castroviejo Blade Breaker
Satin

Catalog #	Overall Length
WI-021-13	13 cm – 5



Catlin Amputation Knife
Satin

Catalog #	Overall Length
WI-021-14	29 cm – 11 1/2



Clappison Carver Fig. 4
Satin

Catalog #	Overall Length
WI-030-07	



Collin Scalpel Handle No. 9
Satin

Catalog #	Overall Length
WI-030-17	13.5 cm – 5 1/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manici per bisturi

Converse Rhinoplastic Knife
 Satin



Catalog #	Overall Length
WI-020-14	16 cm – 6 1/4

Cottle Rhinoplastic Knife
 Satin



Catalog #	Overall Length
WI-021-13	14 cm – 5 1/2

Crane-Kaplan Periodontal Knife Fig. CK-3
 Satin



Fig. CK-3

Catalog #
WI-021-14

Dental Carver
 Satin



Catalog #
WI-030-07


Discoid-Cleoids Carver Fig. 1/2
 Satin



Catalog #
WI-030-17



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Discoid-Cleoids Carver Fig. 3/6
Satin



Catalog #

WI-020-14

Discoid-Cleoids Carver Fig. 4/5
Satin



Catalog #

WI-021-13

Discoid-Cleoids Carver Fig. 89/92
Satin



Catalog #

WI-021-14

Discoid-Cleoids Carver Fig. 90B
Satin



Catalog #

WI-030-07

Discoid-Cleoids Carver Fig. W-A-C-D
Satin



Catalog #

WI-030-17



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Esthetic Carving Knife Fig. 20
 Satin



Catalog #

WI-020-14

Esthetic Carving Knife Fig. 21
 Satin



Catalog #

WI-021-13

Fahnestock Wax Knife Fig. 1
 Mirror



Fig. 1, Wooden Handle

Catalog #

WI-021-14

Fahnestock Wax Knife Fig. 2
 Mirror



Fig. 2, Wooden Handle

Catalog #

WI-030-07

Fisch Footplate Perforator
 Satin



Catalog #

WI-030-17

Overall Length

16 cm – 6 1/4



Scalpel Handles
Skalpelligriffe
Mangos de bisturías
Manches pour bistouris
Manici per bisturi

Fomon Rhinoplastic Knife
Satin



Catalog #	Overall Length
WI-020-14	16.5 cm – 6 1/2

Fowler Tonsil Knife
Satin



Catalog #	Overall Length
WI-021-13	20 cm 8"

Frahm Carver Fig. 2/3
Satin



Catalog #
WI-021-14

Humby – Braithwaite Dermatome (Skin Graft Knife)
Satin




Catalog #	Overall Length
WI-030-07	32 cm – 13

Pudendus Anaesthesia Needle
Satin



Complete With Cannula

Catalog #	Overall Length
WI-030-17	18 cm – 7



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Scalpel Handle No. 7
 Satin



Catalog #	Overall Length
WI-020-14	16 cm – 6 1/4

Silver Dermatome (Skin Graft Knife)
 Satin



Catalog #	Overall Length
WI-021-13	19 cm – 7 1/2

Silver Skin Graft Knife Blade
 Satin



Catalog #
WI-021-14

Virchow Autopsy Knife
 Satin



Catalog #	Overall Length
WI-030-07	25.5 cm – 10

Virchow Cartilage Knife
 Satin



Catalog #	Overall Length
WI-030-17	21cm – 8 1/4



Scalpel Handles
 Skalpgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi



Watson Dermatome (Skin Graft Knife)
 Satin



Catalog #	Overall Length
WI-020-14	30.5 cm – 12

McCoy Laryngoscope
 Satin



Catalog #
WI-021-13

McIntosh Laryngoscope
 Satin




Catalog #
WI-021-14

Miller Laryngoscope
 Satin



Catalog #
WI-030-07

Agarwal Karate Phaco Chopper
 Satin



Catalog #	Overall Length
WI-030-17	11.5 cm – 4 1/2



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris
 Manici per bisturi

Aguja para Cuerpos Extraños Tipo Golf-Club
 Satin

Catalog #	Overall Length
WI-020-14	12.5 cm – 5



Akahoshi Nucleus Sustainer
 Satin

Catalog #	Overall Length
WI-021-13	11.5 cm – 4 1/2



Alfonso Cutting Platform
 Satin

Catalog #	Overall Length
WI-021-14	12.5 cm – 5



Alfonso Nucleus Trisector
 Satin

Catalog #	Overall Length
WI-030-07	10.5 cm – 4 1/4



Alio-Kelman IOL Hook
 Satin

Catalog #	Overall Length
WI-030-17	12 cm – 4 3/4





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Amalgam Carrier
Satin

Catalog #

WI-020-14



Anis Nucleus Manipulator
Satin

Catalog #

WI-021-13

Overall Length

12 cm – 4 3/4



Asa de New Orleans
Satinado, Satin

Catalog #

WI-021-14

Overall Length

13 cm – 5 , 13 cm – 5



Asa de Wilder
Satinado, Satin

Catalog #

WI-030-07

Overall Length

13.5 cm – 5 1/4 , 13.5 cm – 5 1/4



Babinsky Percussion Hammer
Satin

Catalog #

WI-030-17

Overall Length

22 cm – 8 3/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Barraquer Iris Spatula
 Satin

Catalog #	Overall Length
WI-020-14	12 cm – 4 3/4



Bechert Lens Pusher
 Satin

Catalog #	Overall Length
WI-021-13	11.5 cm – 4 1/2



Bechert Nucleus Rotator
 Satin

Catalog #	Overall Length
WI-021-14	11.5 cm – 4 1/2



Binder Flap Elevator
 Satin

Catalog #	Overall Length
WI-030-07	13 cm – 5



Blumenthal Conjunctiva Dissector
 Satin

Catalog #	Overall Length
WI-030-17	13 cm – 5





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Bonn Micro Iris Hook
Satin



Catalog #	Overall Length
WI-020-14	12.5 cm – 5

Brierley Nucleus Splitter
Satin



Catalog #	Overall Length
WI-021-13	12 cm – 4 3/4

Brierley Nucleus Splitter
Satin



Catalog #	Overall Length
WI-021-14	12 cm – 4 3/4

Bucher-Kansas Nucleus Trisector
Satin



Catalog #	Overall Length
WI-030-07	10.5 cm – 4 1/4

Bunge Evisceration Spoon
Satin



Catalog #	Overall Length
WI-030-17	14 cm – 5 1/2



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris
 Manici per bisturi

Buratto Double Ended Lasik Flap Protector
 Satin



Catalog #	Overall Length
WI-020-14	11.5 cm – 4 1/2

Buratto Lasik Flap Protector
 Satin



Catalog #	Overall Length
WI-021-13	12.5 cm – 5

Carter Sphere Introducer
 Satin



Catalog #	Overall Length
WI-021-14	12.5 cm – 5

Castroviejo Cyclodialysis Spatula
 Satin



Catalog #	Overall Length
WI-030-07	13 cm – 5


Castroviejo Synechia Spatula
 Satin



Catalog #	Overall Length
WI-030-17	13 cm – 5



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Chang Micro Finger Chopper
Satin



Catalog #	Overall Length
WI-020-14	11.5 cm – 4 1/2

Chang Phaco Chopper
Satin



Catalog #	Overall Length
WI-021-13	13 cm – 5

Clayman Nucleus Rotator
Satin




Catalog #	Overall Length
WI-021-14	12 cm – 4 3/4

Clayman-Jaffe Lens Manipulating Hook
Satin



Catalog #	Overall Length
WI-030-07	11.5 cm – 4 1/2

Corneal Dissector
Satin



Catalog #	Overall Length
WI-030-17	12 cm – 4 3/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris
 Manici per bisturi

Crawford Retrieval Hook
 Satin



Catalog #	Overall Length
WI-020-14	14.5 cm – 5 3/4

Cucharilla de Cullen
 Satinado, Satin



Catalog #	Overall Length
WI-021-13	13 cm – 5 , 13 cm – 5

Cucharilla de Paton
 Satinado, Satin



Catalog #	Overall Length
WI-021-14	14 cm – 5 1/2 , 14 cm – 5 1/2

Culler Iris Spatula
 Satin



Catalog #	Overall Length
WI-030-07	13 cm – 5

Culler Lens Spoon
 Satin



Catalog #	Overall Length
WI-030-17	13 cm – 5



Scalpel Handles
Skalpelgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Dastoor Lacrimal Sac Dissector
Satin



Catalog #	Overall Length
WI-020-14	13 cm – 5

Dastoor Pupil and Iris Repositor
Satin



Catalog #	Overall Length
WI-021-13	13 cm – 5

Deitz ICL Manipulator
Satin



Catalog #	Overall Length
WI-021-14	12 cm – 4 3/4

Dodick-Kammann Nucleus Chopper
Satin



Catalog #	Overall Length
WI-030-07	12 cm – 4 3/4

Drysdale Nucleus Manipulator
Satin



Catalog #	Overall Length
WI-030-17	11 cm – 4 1/4



Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Ellis Foreign Body Spud
Satin

Catalog #	Overall Length
WI-020-14	11 cm – 4 1/4



Lightly Curved and Rounded Tip

Elschnig Cyclodialysis Spatula
Satin

Catalog #	Overall Length
WI-021-13	13 cm – 5




Escaf Nucleus Divider
Satinado, Satin

Catalog #	Overall Length
WI-021-14	12 cm – 4 3/4




Escaf Nucleus Manipulator
Satin

Catalog #	Overall Length
WI-030-07	12 cm – 4 3/4



Escaf Nucleus Sustainer
Satin

Catalog #	Overall Length
WI-030-17	12 cm – 4 3/4





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Espátula para Iris de Barraquer
Satinado, Satin



Catalog #	Overall Length
WI-020-14	12 cm – 4 3/4 , 12 cm – 4 3/4

Espátula para Sinequia de Girard
Satinado, Satin



Catalog #	Overall Length
WI-021-13	11.5 cm – 4 1/2 , 11.5 cm – 4 1/2

Fenzl Lens Manipulating Hook
Satin




Catalog #	Overall Length
WI-021-14	12 cm – 4 3/4

Fixation Hook
Satin



Catalog #	Overall Length
WI-030-07	13 cm – 5


Microchoper de Chang
Satin



Catalog #	Overall Length
WI-030-17	11.5 cm – 4 1/2 , 11.5 cm – 4 1/2

Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris
 Manici per bisturi

Minardi Phaco Chopper
 Satin



Catalog #	Overall Length
WI-020-14	11.5 cm – 4 1/2

New Orleans Lens Loop
 Satinado, Satin



Catalog #	Overall Length
WI-021-13	13 cm – 5

Otis Urethrotome
 Satin



Catalog #	Overall Length
WI-021-14	

Rosen Phaco Splitter
 Satin



Catalog #	Overall Length
WI-030-07	11.5 cm – 4 1/2

Rowen Combination Lasik Flap Protector and Spatula
 Satin



Catalog #	Overall Length
WI-030-17	13 cm – 5



Scalpel Handles
Skalpellgriffe
Mangos de bisturias
Manches pour bistouris
Manici per bisturi


Siepsner Endocapsular Controller
Satin

Catalog #	Overall Length
WI-020-14	11.5 cm – 4 1/2



Sims Uterine Depressor
Satin

Catalog #	Overall Length
WI-021-13	26 cm – 10 1/4



Sloane LASEK Epi Peeler
Satin

Catalog #	Overall Length
WI-021-14	11.5 cm – 4 1/2



Sloane LASEK Flap Repositor
Satin

Catalog #	Overall Length
WI-030-07	12 cm – 4 3/4



Sloane LASEK Micro Hoe
Satin


Catalog #	Overall Length
WI-030-17	12 cm – 4 3/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Sloane LASEK Micro Hoe & Epi Peeler
 Satin

Catalog #	Overall Length
WI-020-14	13.5 cm – 5 1/4



Snare Wire
 Satin

Catalog #	Overall Length
WI-021-13	10.5 cm – 4



Tennant Nucleus Rotator
 Satinado, Satin

Catalog #	Overall Length
WI-021-14	12 cm – 4 3/4



Trisector de Kansas
 Satinado, Satin

Catalog #	Overall Length
WI-030-07	11 cm – 4 1/4 , 11 cm – 4 1/4



Wilder Lens Loop
 Satin

Catalog #	Overall Length
WI-030-17	13.5 cm – 5 1/4





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi




Adson Needle Holder
Satin

Catalog #	Overall Length
WI-020-14	18 cm – 7



Arruga Micro Needle Holder
Satin

Catalog #	Overall Length
WI-021-13	13.5 cm – 5 1/4



Barraquer Micro Needle Holder
Satin

Catalog #	
WI-021-14	



Barraquer Micro Needle Holder
Satin

Catalog #	
WI-030-07	



Bozeman Needle Holder
Satin

Catalog #	
WI-030-17	



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturías
 Manches pour bistouris
 Manici per bisturi

Bozemann Knot Tier
 Satin



Catalog #	Overall Length
WI-020-14	21 cm – 8 1/4

Brown Needle Holder
 Satin



Catalog #	Overall Length
WI-021-13	

Brunner Ligature Conductor
 Satin



Catalog #	Overall Length
WI-021-14	18 cm – 7

Caspar Ligature Adjuster
 Satin



Catalog #	Overall Length
WI-030-07	23 cm – 9

Castroviejo Micro Needle Holder
 Satin



Catalog #	Overall Length
WI-030-17	



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Chest Tube Passer
Satin

Catalog #	Overall Length
WI-020-14	24 cm – 9 1/2



Childe Approximation Forcep
Satin

Catalog #	Overall Length
WI-021-13	17.5 cm – 7



Claus Van Eicken Ligature Needle
Satin

Catalog #	Overall Length
WI-021-14	25 cm – 9 3/4



Collier Needle Holder
Satin

Catalog #	Overall Length
WI-030-07	12.5 cm – 5



Collin Clip Applying Forcep
Satin

Catalog #	Overall Length
WI-030-17	13 cm – 5



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manici per bisturi

Converse Needle Holder
 Satin

Catalog # _____
 WI-020-14



Cooper Ligature Needle
 Satin

Catalog # _____ Overall Length _____
 WI-021-13 18.5 cm – 7 1/4



Crile Needle Holder
 Satin

Catalog # _____ Overall Length _____
 WI-021-14 15 cm – 6



Crile-Murray Needle Holder
 Satin

Catalog # _____ Overall Length _____
 WI-030-07 15 cm – 6




Crile-Wood Needle Holder
 Satin

Catalog # _____
 WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



De Martel Wire Saw Guide
Satin



Catalog #	Overall Length
WI-020-14	33 cm – 13

DeBakey Needle Holder
Satin



Catalog #
WI-021-13

Demel Wire Guides
Satin



Catalog #
WI-021-14

Derf Needle Holder
Satin



Catalog #	Overall Length
WI-030-07	12 cm – 4 3/4

Deschamps Ligature Needle
Satin



Catalog #
WI-030-17



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi


Finochietto Needle Holder
 Satin

Catalog #	Overall Length
WI-020-14	26.5 cm – 10 1/2



Gillies Needle Holder
 Satin

Catalog #	Overall Length
WI-021-13	16 cm – 6 1/4



Halsey Needle Holder
 Satin

Catalog #
WI-021-14



Heaney Needle Holder
 Satin

Catalog #
WI-030-07



Hegar-Baumgartner Needle Holder
 Satin

Catalog #
WI-030-17





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Hegenbarth Clip Applying Forcep
Satin



Catalog #	Overall Length
WI-020-14	12 cm – 4 3/4

Hoesel Needle Holder
Satin



Catalog #	Overall Length
WI-021-13	25 cm – 9 3/4

Jacobson Micro Needle Holder
Satin




Catalog #	Overall Length
WI-021-14	

Jacobson Micro Suture Pusher
Satin



Catalog #	Overall Length
WI-030-07	18.5 cm – 7 1/4

Johnson Needle Holder
Satin



Catalog #	Overall Length
WI-030-17	27 cm – 10 3/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Jonesco Ligature Needle
 Satin

Catalog #	Overall Length
WI-020-14	18.5 cm – 7 1/4



Kalt Micro Needle Holder
 Satin

Catalog #	Overall Length
WI-021-13	13.5 cm – 5 1/4



Mathieu Needle Holder
 Satin

Catalog #	Overall Length
WI-021-14	



McPherson Micro Needle Holder
 Satin

Catalog #	Overall Length
WI-030-07	10 cm – 4



Reverdin Ligature Needle
 Satin

Catalog #	Overall Length
WI-030-17	





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Alexander Bone Chisel
Dual Satin & Mirror

Catalog #	Overall Length
WI-020-14	17.5 cm – 7



Alexander Bone Gouge
Dual Satin & Mirror

Catalog #	Overall Length
WI-021-13	17.5 cm – 7



Ballenger Septum Gouge
Dual Satin & Mirror

Catalog #	Overall Length
WI-021-14	15 cm – 6



Bauer Osteotome
Dual Satin & Mirror

Catalog #	Overall Length
WI-030-07	20.5 cm – 8



Bone Chisel
Dual Satin & Mirror

Catalog #	Overall Length
WI-030-17	13.5 cm – 5 1/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manici per bisturi

Bone Gouge
 Dual Satin & Mirror

Catalog # Overall Length
 WI-020-14 13.5 cm – 5 1/4



Bone Gouge for Lamina
 Dual Satin & Mirror

Catalog # Overall Length
 WI-021-13 23 cm – 9 1/4



Bone Osteotome
 Dual Satin & Mirror

Catalog # Overall Length
 WI-021-14 13.5 cm – 5 1/4




Brunetti Post Mortem Chisel
 Satin

Catalog # Overall Length
 WI-030-07 28.5 cm – 11 1/4




Buckingham Chisel
 Satin

Catalog # Overall Length
 WI-030-17 18 cm – 7





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Cattoni Hoe/Chisel Fig. 104/105
Satin

Catalog #

WI-020-14



Fig. 104/105

Chisel Fig. 11/12
Satin

Catalog #

WI-021-13



Chisel Fig. 40/41
Satin

Catalog #

WI-021-14



Chisel Fig. 7/10
Satin

Catalog #

WI-030-07



Chisel Fig. 8/9
Satin

Catalog #

WI-030-17



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris
 Manici per bisturi

Cinzel de Cottle

Dual Satin & Mirror, Dual Satin & Mirror

Catalog #	Overall Length
WI-020-14	18.5 cm – 7 1/4 , 18.5 cm – 7 1/4



Cinzel de Lexer

Dual Satin & Mirror, Dual Satin & Mirror

Catalog #	Overall Length
WI-021-13	20 cm – 8 , 20 cm – 8



Cinzel Estándar

Dual Satin & Mirror, Dual Satin & Mirror

Catalog #	Overall Length
WI-021-14	13.5 cm – 5 1/4 , 13.5 cm – 5 1/4



Cinelli Bone Osteotome

Dual Satin & Mirror

Catalog #	Overall Length
WI-030-07	16 cm – 6 1/4



Cloward Bone Gouge

Dual Satin & Mirror

Catalog #	Overall Length
WI-030-17	24.5 cm – 9 3/4





Scalpel Handles
Skalpelligriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi


Cloward Bone Osteotome
Dual Satin & Mirror

Catalog #	Overall Length
WI-020-14	24.5 cm – 9 3/4



Cottle Bone Osteotome
Dual Satin & Mirror

Catalog #	Overall Length
WI-021-13	18.5 cm – 7 1/4



Councilman's Chisel
Satin

Catalog #	Overall Length
WI-021-14	18 cm – 7



Cushing Chisel Scaler Fig. 61
Satin

Catalog #	Overall Length
WI-030-07	



Cushing Chisel Scaler Fig. 62
Satin

Catalog #	Overall Length
WI-030-17	



Scalpel Handles
Skalpelli-griffe
Mangos de bisturíes
Manches pour bistouris

Dahmen Bone Osteotome
Dual Satin & Mirror

Catalog # Overall Length
WI-020-14 30 cm – 12



Darby-Perry Chisel Fig. DP5K/6K
Satinado, Satin

Catalog #
WI-021-13



Fig. DP5K/6K

Darby-Perry Chisel Fig. DP9K/11K
Satin

Catalog #
WI-021-14



Fig. DP9K/11K

Ferrozell Handle Bone Chisel
Dual Satin & Mirror

Catalog #
WI-030-07



Gardner Bone Chisel
Dual Satin & Mirror

Catalog # Overall Length
WI-030-17 17 cm – 6 3/4





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Hibbs Bone Chisel
Dual Satin & Mirror

Catalog #	Overall Length
WI-020-14	24.5 cm – 9 3/4



Hibbs Bone Osteotome
Dual Satin & Mirror

Catalog #	Overall Length
WI-021-13	24.5 cm – 9 3/4




Hoke Bone Osteotome
Dual Satin & Mirror

Catalog #	Overall Length
WI-021-14	17 cm – 6 3/4



Lambotte Bone Osteotome
Dual Satin & Mirror

Catalog #	Overall Length
WI-030-07	24 cm – 9 1/2



Lexer Bone Chisel
Dual Satin & Mirror

Catalog #	Overall Length
WI-030-17	20 cm – 8



Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Lexer Bone Osteotome
Dual Satin & Mirror

Catalog # _____
WI-020-14



Lucas Bone Chisel
Dual Satin & Mirror

Catalog # _____ Overall Length _____
WI-021-13 17 cm – 6 3/4




Lucas Bone Osteotome
Dual Satin & Mirror

Catalog # _____ Overall Length _____
WI-021-14 17 cm – 6 3/4



Mini-Lambotte Bone Osteotome
Dual Satin & Mirror

Catalog # _____
WI-030-07




Partsch Bone Chisel
Dual Satin & Mirror

Catalog # _____
WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Partsch Bone Gouge
Dual Satin & Mirror

Catalog #
WI-020-14



Post Mortem Chisel
Satin

Catalog # Overall Length
WI-021-13 13 cm – 5 1/4



Skull Breaker & Chisel
Satin

Catalog # Overall Length
WI-021-14 14 cm – 5 1/2



Smith-Peterson Bone Osteotome
Dual Satin & Mirror

Catalog # Overall Length
WI-030-07 20.5 cm – 8



Snyder OSU Bone Gouge
Dual Satin & Mirror

Catalog # Overall Length
WI-030-17 16.5 cm – 6 1/2



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Stille Bone Chisel
 Dual Satin & Mirror

Catalog #	Overall Length
WI-020-14	20 cm – 8




Stille Bone Osteotome
 Dual Satin & Mirror

Catalog #	Overall Length
WI-021-13	20 cm – 8



Wagner Bone Gouge
 Dual Satin & Mirror

Catalog #	Overall Length
WI-021-14	34 cm – 13 1/2



Walter Bone Chisel
 Dual Satin & Mirror

Catalog #	Overall Length
WI-030-07	19 cm – 7 1/2



Walter Bone Osteotome
 Dual Satin & Mirror

Catalog #	Overall Length
WI-030-17	19 cm – 7 1/2





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Adams Utility Plier
Satin

Catalog # _____
WI-020-14



Adhesive Removing Plier
Satin

Catalog # _____
WI-021-13



Arch Bending Plier
Satin

Catalog # _____
WI-021-14



Band Removing Plier
Satin

Catalog # _____
WI-030-07



Band Seating Plier
Satin

Catalog # _____
WI-030-17



Scalpel Handles
Skalpellgriffe
Mangos de bisturís
Manches pour bistouris
Manici per bisturi

Bird Beak Plier
Satin

Catalog #

WI-020-14



Bird Beak Plier With Cutter
Satin

Catalog #

WI-021-13



Bracket Removing Plier
Satinado, Satin

Catalog #

WI-021-14



Combi-Flat Nose Gripping & Cutting Pliers
Satin

Catalog #

WI-030-07



Coon Suture Tying Plier
Satin

Catalog #

WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi




Corwin Wire Twisting Forceps
Satin

Catalog #	Overall Length
WI-020-14	16 cm – 6 1/4



Crown & Band Contouring Plier
Satin

Catalog #
WI-021-13



De La Rosa Hollow Chop Contouring Plier
Satin

Catalog #
WI-021-14



Demel Wire Tightener and Twister
Satin

Catalog #	Overall Length
WI-030-07	27 cm – 10 3/4



Flush Cut & Hold Distal End Cutter
Satin

Catalog #
WI-030-17



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris
 Manici per bisturi

Flush Cut Distal End Cutter
 Satin

Catalog # _____
 WI-020-14



Force Module Separating Plier
 Satin

Catalog # _____
 WI-021-13



Front Cutting Action Pin and Wire Cutters
 Satin

Catalog # _____
 WI-021-14



Hammerhead Tie Back Plier
 Satin

Catalog # _____
 WI-030-07



Heavy Duty Plate, Pin and Wire Cutters
 Satin

Catalog # _____
 WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Hollow Chop Plier
Satin

Catalog # _____
WI-020-14



Hook Crimping Plier
Satin

Catalog # _____
WI-021-13



How Utility Plier
Satin

Catalog # _____
WI-021-14



Straight, Curved

Jaraback Wire Forming Plier
Satin

Catalog # _____
WI-030-07



Lateral Cutting Action Pin and Wire Cutters
Satin

Catalog # _____ Overall Length
WI-030-17 23 cm – 9



For Soft Wires Up to 3.0 mm Ø, For Hard Wires Up to 2.2 mm Ø



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici de bisturi

Ligature Cutter
Satin

Catalog #
WI-020-14



Light Wire Bird Beak Plier
Satin

Catalog #
WI-021-13



Lingual Arch Bending Plier
Satin

Catalog #
WI-021-14



Long Tapered Bird Beak Plier
Satin

Catalog #
WI-030-07



Loop Forming Plier
Satin

Catalog #
WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Loute Wire Tightener
Satin



Catalog # Overall Length
WI-020-14 21 cm – 8 1/4

Marcotte Loop Forming Plier
Satin



Catalog #
WI-021-13

Micro Mini Wire & Ligature Cutter
Satin



Catalog #
WI-021-14

Mini Three Jaw Plier
Satin



Catalog #
WI-030-07

Mini Wire & Ligature Cutter
Satin



Catalog #
WI-030-17



Scalpel Handles
 Skalpelli-griffe
 Mangos de bisturís
 Manches pour bistouris
 Manici ner bisturi

Multi Use Wire & Ligature Cutter
 Satin



Catalog #

WI-020-14

Nance Loop Forming Plier
 Satin



Catalog #

WI-021-13


O'Brien Orthodontic Plier
 Satin



Catalog #

WI-021-14

Obwegeser Wire Loop Forceps
 Satin



Catalog #

WI-030-07

Overall Length

15 cm – 6

Omega Loop Forming Plier
 Satin



Catalog #

WI-030-17



Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Optical Plier
Satin

Catalog # _____
WI-020-14



Orthodontic Finishing Step Plier
Satin

Catalog # _____
WI-021-13



Parallel Grip Flat Nose Plier
Satin

Catalog # _____ Overall Length
WI-021-14 18.5 cm – 7 1/4



Pin and Wire Cutters
Satin

Catalog # _____
WI-030-07



Side Wire Cutters
Satin

Catalog # _____ Overall Length
WI-030-17 13 cm – 5 1/4



For Soft Wires Up to 1.5 mm Ø. For Hard Wires Up to 1.0 mm Ø



Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manis para bisturí

Three Jaw Plier
Satin



Catalog #

WI-020-14

Tip Back Plier
Satin



Catalog #

WI-021-13


Tweed Loop Forming Plier
Satin



Catalog #

WI-021-14

Wire Twisting Forceps
Satin



Catalog #

WI-030-07

Young Loop Forming Plier
Satin



Catalog #


WI-030-17



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Bakes Gall Duct Dilators
Satin

Catalog #	Overall Length
WI-020-14	32 cm – 12 1/2



Bowman Probes
Satin

Catalog #	Overall Length
WI-021-13	13 cm – 5



Bridge Artery Forcep
Satin

Catalog #	Overall Length
WI-021-14	27.5 cm – 10 3/4



Brodie Butterfly Probe / Grooved Director
Satin

Catalog #	Overall Length
WI-030-07	20 cm – 8



Butterfly Probe / Grooved Director
Satin

Catalog #	Overall Length
WI-030-17	



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturías
 Manches pour bistouris
 Manici per bisturi

Castroviejo Lacrimal Dilator
 Satin



Catalog #	Overall Length
WI-020-14	14 cm – 5 1/2

Cawthorne Probe
 Satin



Catalog #	Overall Length
WI-021-13	14 cm – 5 1/2

Cooley Vascular Dilator
 Satin



Malleable

Catalog #	Overall Length
WI-021-14	13 cm – 5

Cotton Applicator
 Satin



Round Spiral End End

Catalog #	Overall Length
WI-030-07	18 cm – 7

Cowhorn Dental Explorer Fig. 3CH
 Satin




Fig. 3CH Composite™

Fig. 3CH

Catalog #
WI-030-17




Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Crawford Lacrimal Intubation Set
Satin

Catalog # _____
WI-020-14



DeBakey Vascular Dilator
Satin

Catalog # _____ Overall Length _____
WI-021-13 19 cm – 7 1/2



Malleable

Dental Explorer Fig. 2R/2L
Satin

Catalog # _____
WI-021-14



Fig. 2R/2L

Dental Explorer Fig. 3
Satin

Catalog # _____
WI-030-07



Fig. 2R/2L

Dental Explorer Fig. 3A
Satin

Catalog # _____
WI-030-17



Fig. 3A



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Double Ended Dental Explorer Fig. 3
 Satin



Catalog # _____
 WI-020-14

Double Ended Probe
 Satin



Catalog # _____
 WI-021-13

Doyen Butterfly Probe / Grooved Director
 Satin



Catalog #	Overall Length
WI-021-14	14.5 cm – 5 3/4

Explorer Fig. Dg16
 Satin



Catalog # _____
 WI-030-07

Explorer Fig. DG16/17
 Satin



Catalog # _____
 WI-030-17



Scalpel Handles
Skalpelgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi



Explorer Fig. DG16/23
Satin

Catalog # _____
WI-020-14



Fig. DG16/23

Farrell Cotton Applicator
Satin

Catalog # _____
WI-021-13



Triangular Shaped End
Round Spiral End

Fillion Explorer
Satin

Catalog # _____
WI-021-14



Garret Vascular Dilator
Satin

Catalog # _____
WI-030-07



Malleable

Glickman Measuring Probe Fig. 26 G
Satin


Catalog # _____
WI-030-17



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Goldman-Fox Explorer/Measuring Probe Fig. 23/GF
 Satin

Catalog # _____
 WI-020-14




Goldman-Fox Measuring Probe
 Satin

Catalog # _____
 WI-021-13



Goldman-Fox Williams Measuring Probe
 Satin

Catalog # _____
 WI-021-14



Goodell Uterine Dilator
 Satin

Catalog # Overall Length
 WI-030-07 34 cm – 13 1/2



With graduations and fixing screw

Hank Uterine Dilators
 Satin

Catalog # Overall Length
 WI-030-17 26.5 cm – 10 1/2



Set of 6



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Harms Trabeculotomy Probe
Satin

Catalog #	Overall Length
WI-020-14	5 cm – 2



Hartmann Cotton Applicator
Satin

Catalog #	Overall Length
WI-021-13	22 cm – 8 1/2



Hegar Uterine Dilators
Mirror

Catalog #	Overall Length
WI-021-14	




Infant Lacrimal Dilator
Satin

Catalog #	Overall Length
WI-030-07	8 cm – 3



Jacobson Micro Vessel Probe
Satin

Catalog #	Overall Length
WI-030-17	18.5 cm – 7 1/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris
 Manici per bisturi

Jobson-Horne Ear Probe
 Satin

Catalog # _____
 WI-020-14



Killian Probe
 Satin

Catalog # _____ Overall Length _____
 WI-021-13 16.5 cm – 6 1/2



Kuhn Frontal Ostium Seeker
 Satin

Catalog # _____ Overall Length _____
 WI-021-14 22 cm – 8 1/2




Laborde Tracheal Dilator
 Satin

Catalog # _____ Overall Length _____
 WI-030-07 14 cm – 5 1/2



Lockhart-Mummery Fistula Probe
 Satin

Catalog # _____ Overall Length _____
 WI-030-17 16.5 cm – 6 1/2






Scalpel Handles
Skalpelligriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Martin Vascular Tunneling Set
Satin

Catalog # _____

WI-020-14

- Set Contents**
- 1 - Handle for Tunneling Set
 - 2 - Tunneling Shafts x 3
 - 3 - Tunneling Tip x 4 (6 mm, 8 mm, 10 mm, 12 mm)
 - 4 - Vein Tip x 1
 - 5 - Luer Adapter
 - 6 - Storage Case Stainless



Maxillary Ostium Sinus Seeker
Satin

Catalog # _____

WI-021-13

Overall Length

19 cm – 7 1/2



Mayo Common Duct Probe
Satin

Catalog # _____

WI-021-14

Overall Length

25 cm – 10



Measuring Probe Fig. CP-10
Satin

Catalog # _____

WI-030-07



Measuring Probe Fig. CP-11
Satin

Catalog # _____

WI-030-17



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 sturi

Measuring Probe Fig. CP-12
 Satin

Catalog # _____
 WI-020-14



Measuring Probe Fig. CP-2
 Satin

Catalog # _____
 WI-021-13



Measuring Probe Fig. CP-8
 Satin

Catalog # _____
 WI-021-14



Nail Probe
 Satin

Catalog #	Overall Length
WI-030-07	14 cm – 5 1/2



Nail Spatula and Probe
 Satin

Catalog #	Overall Length
WI-030-17	16 cm – 6 1/4





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Abdominal Spatula
Satin



Catalog #	Overall Length
WI-020-14	33 cm – 13

Adson (Cone) Self Retaining Retractor
Satin



Catalog #	Overall Length
WI-021-13	

Adson Dura Hook
Satin



Catalog #	Overall Length
WI-021-14	20 cm – 8

Adson Hinged Self Retaining Retractor
Satin



Catalog #	Overall Length
WI-030-07	33 cm – 13

Adson Self Retaining Retractor
Satin



Catalog #	Overall Length
WI-030-17	26 cm – 10 1/4 , 29 cm – 11 1/2



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manches per bisturi

Agricola Lacrimal Sac Retractor
 Satin

Catalog # _____
 WI-020-14



Alio-Kelman IOL Hook
 Satin

Catalog # _____ Overall Length _____
 WI-021-13 12 cm – 4 3/4



Allison Lung Spatula
 Satin

Catalog # _____
 WI-021-14



Allport Wound Spreader
 Satin

Catalog # _____ Overall Length _____
 WI-030-07 9.5 cm – 3 3/4



Alm Wound Spreader
 Satin

Catalog # _____
 WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi


Anderson Claw Retractor
Satin

Catalog #	Overall Length
WI-020-14	12 cm – 4 3/4



Anderson-Adson Self Retaining Retractor
Satin

Catalog #	Overall Length
WI-021-13	20 cm – 8



Balfour-Baby Self Retaining Retractor
Satin

Catalog #	Overall Length
WI-021-14	



Barraquer Iris Spatula
Satin

Catalog #	Overall Length
WI-030-07	12 cm – 4 3/4



Beeson Plaster Cast Spreader
Satin

Catalog #	Overall Length
WI-030-17	30.5 cm – 12



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manici per bisturi

Binder Flap Elevator
 Satin

Catalog #	Overall Length
WI-020-14	13 cm – 5




Blount Knee Retractor
 Satin

Catalog #	Overall Length
WI-021-13	17.5 cm – 7



Blount Retractor
 Satin

Catalog #	Overall Length
WI-021-14	26.5 cm – 10 1/2



Bonn Micro Iris Hook
 Satin

Catalog #	Overall Length
WI-030-07	12.5 cm – 5



Bose Fat Spreader
 Satin

Catalog #	Overall Length
WI-030-17	6 cm – 2 1/4





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Bose Tracheal Retractor
Satin

Catalog #	Overall Length
WI-020-14	16 cm – 6 1/4




Bose Tracheal Spreader
Satin

Catalog #	Overall Length
WI-021-13	6 cm – 2 1/4



Bowdler-Henry Retractor
Satin

Catalog #	Overall Length
WI-021-14	17 cm – 6 3/4



Brain Spatula
Satin

Catalog #	Overall Length
WI-030-07	20 cm – 8



Malleable
12 mm / 17 mm

Brom Vein Hook
Satin

Catalog #	Overall Length
WI-030-17	19.5 cm – 7 3/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris
 Manici per bisturi

Bruenings Tongue Depressor
 Satin

Catalog #	Overall Length
WI-020-14	19 cm – 7 1/2



Brunner Retractor
 Satin

Catalog #	Overall Length
WI-021-13	25.5 cm – 10




Buchs Retractor
 Dual Satin & Mirror

Catalog #	Overall Length
WI-021-14	20.5 cm – 8



Buckley Wound Spreader
 Satin

Catalog #	Overall Length
WI-030-07	6 cm – 2 1/4



Bugge Takedown I.M.A Mammaria Retractor
 Dual Satin & Mirror

Catalog #
WI-030-17





Scalpel Handles
Skalpelligriffe
Mangos de bisturías
Manches pour bistouris
Manici per bisturi

Burford Rib Spreader
Dual Satin & Mirror

Catalog # _____
WI-020-14



With 2 Pairs of Lateral Blades, Aluminium

Campbell Nerve Root Retractor
Satin

Catalog # _____ Overall Length _____
WI-021-13 21.5 cm – 8 1/2



Castroviejo Cyclodialysis Spatula
Satin

Catalog # _____ Overall Length _____
WI-021-14 13 cm – 5



Castroviejo Synechiaie Spatula
Satin

Catalog # _____ Overall Length _____
WI-030-07 13 cm – 5



Clayman-Jaffe Lens Manipulating Hook
Satin

Catalog # _____ Overall Length _____
WI-030-17 11.5 cm – 4 1/2



Scalpel Handles
 Skalpgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manici per bisturi

Crawford Retrieval Hook
 Satin



Catalog #	Overall Length
WI-020-14	14.5 cm – 5 3/4

Culler Iris Spatula
 Satin



Catalog #	Overall Length
WI-021-13	13 cm – 5

Elschnig Cyclodialysis Spatula
 Satin



Catalog #	Overall Length
WI-021-14	13 cm – 5

Espátula para Iris de Barraquer
 Satin



Catalog #	Overall Length
WI-030-07	12 cm – 4 3/4 , 12 cm – 4 3/4

Espátula para Sinequia de Girard
 Satinado, Satin



Catalog #	Overall Length
WI-030-17	11.5 cm – 4 1/2 , 11.5 cm – 4 1/2



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Fenzl Lens Manipulating Hook
Satin

Catalog #	Overall Length
WI-020-14	12 cm – 4 3/4



Fixation Hook
Satin

Catalog #	Overall Length
WI-021-13	13 cm – 5



Hennig Plaster Cast Spreader
Satin

Catalog #	Overall Length
WI-021-14	28 cm – 11




Kallmorgen Vaginal Speculum and Retractor
Satin

Catalog #	
WI-030-07	



Kartz-Barraquer Wire Speculum
Satin

Catalog #	
WI-030-17	



Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Kaye Facelift Hook
Satin

Catalog #	Overall Length
WI-020-14	13 cm – 5



Kocher-Langenbeck Retractor
Dual Satin & Mirror

Catalog #	Overall Length
WI-021-13	21.5 cm – 8 1/2



Kristeller Vaginal Speculum and Retractor
Satin

Catalog #	Overall Length
WI-021-14	



Nail Spatula
Satin

Catalog #	Overall Length
WI-030-07	11.5 cm – 4 1/2



Quattro Plaster Cast Spreader
Satin

Catalog #	Overall Length
WI-030-17	23 cm – 9





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Adson Cranial Rongeur
Satin

Catalog #	Overall Length
WI-020-14	20 cm – 8



Ainsworth Rubber Dam Punch
Satin

Catalog #	Overall Length
WI-021-13	16 cm – 6 3/4



Alexander Uterine Biopsy Punch
Satin

Catalog #	Shaft Length
WI-021-14	19.5 cm – 7 3/4



Bacon Cranial Rongeur
Satin

Catalog #	Overall Length
WI-030-07	21 cm – 8 1/4



Baggish Uterine Biopsy Punch
Satin

Catalog #	Shaft Length
WI-030-17	22 cm – 8 3/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manici per bisturi

Bethune Rib Shear
 Satin

Catalog # _____
 WI-020-14



Beyer Bone Rongeur
 Satin

Catalog # _____ Overall Length _____
 WI-021-13 18 cm – 7




Beyer Mastoid Rongeur
 Satin

Catalog # _____
 WI-021-14




Blakesley Nasal Punch
 Satin

Catalog # _____ Shaft Length _____
 WI-030-07 12 cm – 4 3/4



Blumenthal Bone Rongeur
 Satin

Catalog # _____
 WI-030-17





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Bohler Bone Rongeur
Satin

Catalog #	Overall Length
WI-020-14	15.5 cm – 6



Bone Shear
Satin

Catalog #	Overall Length
WI-021-13	23 cm – 9



Brunner Rib Shear
Satin

Catalog #	Overall Length
WI-021-14	32 cm – 12 1/2



Burke Uterine Biopsy Punch
Satin

Catalog #	Shaft Length
WI-030-07	22 cm – 8 3/4



Camellin Lasek Trepine
Satin

Catalog #	Overall Length
WI-030-17	12 cm – 4 3/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturías
 Manches pour bistouris
 Manici per bisturi

Camellin Lasek Trepine & Alcohol Retaining Well
 Satin

Catalog #	Overall Length
WI-020-14	12.5 cm – 5



Chevalier-Jackson Basket Punch
 Satin

Catalog #
WI-021-13



Chevalier-Jackson Square Punch
 Satin

Catalog #
WI-021-14



Christensen Punctum Punch
 Satin

Catalog #	Overall Length
WI-030-07	10 cm – 4



Citelli Mastoid Rongeur
 Satin

Catalog #	Overall Length
WI-030-17	17 cm – 6 3/4






Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Cleveland Bone Rongeur
Satin

Catalog # _____
WI-020-14



Cone Cranial Punch
Satin

Catalog # _____ Overall Length _____
WI-021-13 22.5 cm – 8 3/4



Crescent Jaw Nail Cutter
Satin

Catalog # _____ Overall Length _____
WI-021-14 16 cm – 6 1/4



Cushing Laminectomy Rongeur
Satin

Catalog # _____
WI-030-07



Cuticle Cutter
Satin

Catalog # _____ Overall Length _____
WI-030-17 10 cm – 4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manicini per bisturi

Dahlgren Cranial Rongeur
 Satin

Catalog #	Overall Length
WI-020-14	20 cm – 8



De Vilbiss Cranial Rongeur
 Satin

Catalog #	Overall Length
WI-021-13	20.5 cm – 8



Dermal Punch
 Satin

Catalog #	Overall Length
WI-021-14	5.5 cm – 2 1/4 , 8 cm – 3 1/4



Elliot Scleral Trephine
 Satin

Catalog #	Overall Length
WI-030-07	3.5 cm – 1 1/2



Elliot Trephine Handle
 Satin

Catalog #	Overall Length
WI-030-17	4.5 cm – 1 3/4





Scalpel Handles
Skalpellgriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Endolaryngeal Punch
Satin

Catalog # _____
WI-020-14



Faure Uterine Biopsy Punch
Satin

Catalog # _____ Overall Length _____
WI-021-13 21 cm – 8 1/4



Ferris-Smith Laminectomy Rongeur
Satin

Catalog # _____ Shaft Length _____
WI-021-14 18 cm – 7




Fine Jaw Nail Cutter
Satin

Catalog # _____ Overall Length _____
WI-030-07 13 cm – 5



Foraminotomy Kerrison Punch
Satin


Catalog # _____ Shaft Length _____
WI-030-17 20 cm – 8



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 uri

Friedmann Bone Rongeur
 Satin

Catalog # _____
 WI-020-14



Frykholm Bone Rongeur
 Satin

Catalog # _____ Overall Length _____
 WI-021-13 24.5 cm – 9 3/4



Gass Corneoscleral Punch
 Satin

Catalog # _____
 WI-021-14



Gellhorn Uterine Biopsy Punch
 Satin

Catalog # _____ Overall Length _____
 WI-030-07 24 cm – 9 1/2



Giertz-Stille Rib Shear
 Satin

Catalog # _____ Overall Length _____
 WI-030-17 24.5 cm – 9 3/4





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Gluck Rib Shear
Satin

Catalog #	Overall Length
WI-020-14	21 cm – 8 1/4



Drawing Cut

Hajek Sphenoid Punch
Satin

Catalog #	Overall Length
WI-021-13	14 cm – 5 1/2



Hajek-Kofler Sphenoid Punch
Satin

Catalog #	Overall Length
WI-021-14	14 cm – 5 1/2



Heavy Duty Nail Cutter
Satin

Catalog #	Overall Length
WI-030-07	14 cm – 5 1/2




Ingrown Nail Cutter
Satin

Catalog #	Overall Length
WI-030-17	14.5 cm – 5 3/4




Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manici per bisturi

Keyes Dermal Punch
 Satin



Catalog #	Overall Length
WI-020-14	11 cm – 4 1/2

Nail Cutter
 Satin



Catalog #	
WI-021-13	

Narrow Blade Nail Cutter
 Satin



Catalog #	Overall Length
WI-021-14	13 cm – 5

Stille Plaster of Paris Shear
 Satin



Catalog #	
WI-030-07	

Wire Spring Nail Cutter
 Satin



Catalog #	Overall Length
WI-030-17	13 cm – 5




Scissors

Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Adventitia Micro Scissor Satin




Catalog #	Overall Length
WI-010-03	14.5 cm – 5 3/4

Barraquer Iris Scissor Satin



Catalog #	Overall Length
WI-010-04	5.5 cm – 2 1/4

Barre Vascular Scissor Satin




Catalog #	Size
WI-010-34	28 cm – 11

Beckmann Nasal Scissor Satin



Catalog #	Size
WI-012-03	16.5 cm – 6 1/2

BeeBee Crown Scissor Satin



Catalog #	
WI-020-13	



Scalpel Handles
Skalpgriffe
Manos de bisturís
pour bistouris
per bisturi

Bellucci Micro Scissors
Satin

Catalog #	Size
WI-020-14	8 cm – 3



Bergmann Bandage Scissor
Satin

Catalog #	Size
WI-021-13	23.5 cm – 9 1/4



Biemer Vein Scissor
Satin

Catalog #	Size
WI-021-14	13 cm – 5



Bone Shear
Satin

Catalog #	Size
WI-030-07	23 cm – 9



Bonn Iris Scissor
Satin

Catalog #	Size
WI-030-17	9 cm – 3 1/2





Scalpel Handles
Skalpelli
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi


Bowel Dissecting Scissor
Satin

Catalog #	Size
WI-010-03	28 cm – 11



Bozemann Gynecological Scissor
Satin

Catalog #	Size
WI-010-04	22 cm – 8 3/4



Braun-Stadler Episiotomy Scissor
Satin

Catalog #	Size
WI-010-34	



Bruns Plaster of Paris Shear
Satin

Catalog #	Size
WI-012-03	24 cm – 9 1/2



Buck Stitch Scissor
Satin

Catalog #	Size
WI-020-13	



Scalpel Handles
Skalpelligriffe
Mangos de bisturios
Manches pour bistouris
Manici per bisturi

Busch Umbilical Cord Scissor
Satin

Catalog #	Size
WI-020-14	16 cm – 6 1/4



Caplan Septum Scissor
Satin

Catalog #	Size
WI-021-13	20 cm – 8



Castroviejo Corneoscleral Scissor
Satin

Catalog #	Size
WI-021-14	



Castroviejo Keratoplasty Scissor
Satin

Catalog #	Size
WI-030-07	11.5 cm – 4 1/2



Castroviejo Universal Corneal Scissor
Satin

Catalog #	Size
WI-030-17	11 cm – 4 1/2





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Chadwick Gum Scissor
Satin

Catalog #	Size
WI-010-03	11.5 cm – 4 1/2




Cinelli Dissecting Scissor
Satin

Catalog #	Size
WI-010-04	11.5 cm – 4 1/2




Converse (Graefe) Nasal Scissor
Satin

Catalog #	Size
WI-010-34	10.5 cm – 4 1/4




Cooley Vascular Scissor
Satin

Catalog #	Size
WI-012-03	19 cm – 7 1/2



Cottle Bulldog Scissor
Satin


Catalog #	Size
WI-020-13	11.5 cm – 4 1/2



Scalpel Handles
Skalpelligriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi


Cottle Nasal Scissor
Satin

Catalog #	Size
WI-020-14	16 cm – 6 1/4



Cottle-Masing Plastic Surgery Scissor
Satin

Catalog #	Size
WI-021-13	10.5 cm – 4 1/4



Dean Tonsil Scissor
Satin

Catalog #	Size
WI-021-14	17.5 cm – 7



Deaver Operating Scissor
Satin

Catalog #	Size
WI-030-07	14 cm – 5 1/2



Debakey Arteriotomy Scissor
Satin

Catalog #	Size
WI-030-17	17.5 cm – 7





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

DeBakey Offset Scissor
Satin

Catalog #	Size
WI-010-03	15.5 cm – 6



Delicate Nail Scissor
Satin

Catalog #	Size
WI-010-04	8.5 cm – 3 1/4 , 9 cm – 3 1/2




Diethrich Vascular Scissor
Satin

Catalog #	Size
WI-010-34	17.5 cm – 7



Diethrich-Potts Vascular Scissor
Satin

Catalog #	Size
WI-012-03	18 cm – 7



Enterotomy Scissor
Satin

Catalog #	Size
WI-020-13	21.5 cm – 8 1/2



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturios
 Manches pour bistouris

Esmarch Bandage Scissor
 Satin

Catalog # _____
 WI-020-14



Lister Bandage Scissor
 Satin

Catalog # _____
 WI-021-13



Lister Excentric Bandage Scissor
 Satin

Catalog # _____
 WI-021-14



Lloyd-Davie (Goligher) Rectal Scissor
 Satin

Catalog #	Size
WI-030-07	27 cm – 10 3/4



Lorenz Bandage Scissor
 Satin

Catalog #	Size
WI-030-17	23 cm – 9





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Mayo Dissecting Scissor
Satin

Catalog # _____

WI-010-03




Stainless Steel Blade

Tungsten Carbide Blade

Mayo SuperEdge™ Dissecting Scissor
Satin

Catalog # _____

WI-010-04



Straight Beveled Blades

Curved Beveled Blades




Mueller Rectal Scissor
Satin

Catalog # _____

WI-010-34

Size


32.5 cm – 12 3/4



Nail Scissor
Satin

Catalog # _____

WI-012-03



Straight

Curved



Rectal Scissor
Satin

Catalog # _____

WI-020-13

Size

32.5 cm – 12 3/4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manici per bisturi


Seutin Plaster of Paris Shear
 Satin

Catalog #	Size
WI-020-14	23 cm – 9



Stelzner Rectal Scissor
 Satin

Catalog #	Size
WI-021-13	31 cm – 12 1/4



Systrunk Dissecting Scissor
 Satin

Catalog #	Size
WI-021-14	13 cm – 5



Ultra Delicate Nail Scissor
 Satin

Catalog #	Size
WI-030-07	9 cm – 3 1/2



Universal Bandage Scissor
 Satin

Catalog #
WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Alan-Parks Rectal Speculum
Satin

Catalog #	Overall Length
WI-010-03	14 cm – 5 1/2



With Fiber Optic Illumination

Alphonso Infant Lid Speculum
Satin

Catalog #	
WI-010-04	



With Finger Grip

Auvard Vaginal Speculum
Satin

Catalog #	Size
WI-010-34	23.5 cm – 9 1/4



With Detachable Weight


Barr Rectal Speculum
Satin

Catalog #	Overall Length
WI-012-03	17 cm – 6 3/4



Barraquer Infant Wire Speculum
Satin

Catalog #	
WI-020-13	




9 mm



Scalpel Handles
Skalpgriffe
Mangos de bisturíes
lanches pour bistouris
lanici per bisturi

Barraquer Wire Speculum
Satin

Catalog # _____
WI-020-14



Breisky Vaginal Speculum
Satin

Catalog # _____
WI-021-13



Castroviejo Eye Speculum Adjustable
Satin

Catalog # _____
WI-021-14




Collin Vaginal Speculum
Satin

Catalog # _____
WI-030-07



Cook Eye Speculum
Satin

Catalog # _____
WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Cottle Nasal Speculum
Satin

Catalog #	Overall Length
WI-010-03	14 cm – 5 1/2



Without Screw Fixation and With Screw Fixation

Cushing-Landolt Speculum
Satin

Catalog #	Overall Length
WI-010-04	



75 x 18 mm / 82 x 19 mm / 110 x 19 mm

DeLee Universal Bladder Retractor
Satin

Catalog #	Overall Length
WI-010-34	24 cm – 9 1/2



Doyen Vaginal Speculum
Satin

Catalog #	Overall Length
WI-012-03	



Slightly Concave

Farrion Ear Speculum
Satin

Catalog #	Overall Length
WI-020-13	3.6 cm – 1 1/2



Black



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturís
 Manches pour bistouris
 Manici per bisturi




20 mm
 Solid Blades

Feaster Wire Speculum
 Satin


Catalog # _____
 WI-020-14

Flieringa Scleral Fixation Ring
 Satin

Catalog # _____
 WI-021-13




Set of 9 Ring



Frankel Infant Lid Speculum
 Satin

Catalog # _____
 WI-021-14



Grave Vaginal Speculum
 Satin

Catalog # _____
 WI-030-07




Hartmann-Halle Nasal Speculum
 Satin

Catalog # _____
 WI-030-17






Scalpel Handles
Skalpellsgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Hartmann-Halle Nasal Speculum
Satin

Catalog #	Overall Length
WI-010-03	3.6 cm – 1 1/2



Hirschman Anoscope
Satin

Catalog #
WI-010-04



Infant Lid Speculum
Satin

Catalog #
WI-010-34



Ives-Fansler Anoscope
Satin

Catalog #	Shaft Length
WI-012-03	95 mm



Jackson Vaginal Speculum
Satin

Catalog #
WI-020-13



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi


Kallmorgen Vaginal Speculum and Retractor
 Satin

Catalog # _____
 WI-020-14



Kartz-Barraquer Wire Speculum
 Satin

Catalog # _____
 WI-021-13



Kelly Proctoscope
 Satin

Catalog # _____
 WI-021-14



Kelly Sphincteroscope
 Satin

Catalog # _____
 WI-030-07



Killian Nasal Speculum
 Satin

Catalog # _____ Overall Length
 WI-030-17 13 cm – 5





Scalpel Handles
Skalpgriffe
Mangos de bisturís
Manches pour bistouris
Manici per bisturi

Kogan Endocervical Speculum
Satin

Catalog #	Overall Length
WI-010-03	24.5 cm – 9 3/4



Kristeller Vaginal Speculum and Retractor
Satin

Catalog #
WI-010-04



Lancaster Eye Speculum
Satin

Catalog #
WI-010-34



Landau Vaginal Speculum
Satin

Catalog #	Overall Length
WI-012-03	29 cm – 11 1/2



Landolt Spreader
Satin

Catalog #	Overall Length
WI-020-13	21.5 cm – 8 1/2



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi



Landon Vaginal Speculum
 Satin

Catalog #	Overall Length
WI-020-14	19.5 cm – 7 1/2

Lester-Burch Eye Speculum
 Satin

Catalog #
WI-021-13



Lieberman Eye Speculum
 Satin

Catalog #
WI-021-14



Martin Nasal Speculum
 Satin

Catalog	Overall Length
WI-030-07	12.5 cm – 5 , 13 cm – 5




Mathieu Rectal Speculum
 Satin

Catalog #	Overall Length
WI-030-17	19.5 cm – 7 3/4





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi



Small I.D. 15 mm

Medium I.D. 16 mm

Large I.D. 17 mm

McNeill-Goldman Scleral Fixation Ring & Blepharostat
Satin

Catalog # _____

WI-010-03

Mellinger Eye Speculum
Satin



Catalog # _____

WI-010-04


Mod. Vienna Nasal Speculum
Satin



Catalog # _____

WI-010-34

Murdoch Eye Speculum
Satin



12 mm

15 mm

Catalog # _____

WI-012-03


Parks Rectal Speculum
Satin



With Fiber Optic Illumination

Catalog # _____

WI-020-13



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 ur bistouris
 isturi

Pederson Vaginal Speculum
 Satin

Catalog # _____
 WI-020-14



Pratt Rectal Speculum
 Satin

Catalog # Overall Length
 WI-021-13 20.5 cm – 8



Richnau-Holmgren Ear Speculum
 Satin

Catalog # _____
 WI-021-14



Semm Vaginal Speculum
 Satin

Catalog # _____
 WI-030-07



Siegle Otoscope with Ear Speculums
 Satin

Catalog # _____
 WI-030-17





Scalpel Handles
Skalpellsgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi


Adson Suction Tube
Satin

Catalog #	Overall Length
WI-010-03	20 cm – 8



Air Injection Cannula
Mirror

Catalog #	Overall Length
WI-010-04	



Alio LASIK Irrigation Cannula
Mirror

Catalog #	Overall Length
WI-010-34	



Andrews-Pynchon Suction Tube
Satin

Catalog #	Overall Length
WI-012-03	23 cm – 9



Aspiration Handpiece
Satin

Catalog #	Overall Length
WI-020-13	10.5 cm – 4



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi


Atkinson Retrobulbar Needle
 Mirror

Catalog # _____
 WI-020-14



Banaji LASIK Cannula
 Mirror

Catalog # _____
 WI-021-13



Baron-Schuknecht Suction Tube
 Satin

Catalog #	Shaft Length
WI-021-14	75 mm



Bishop-Harmon Anterior Chamber Irrigator
 Mirror


Catalog # _____
 WI-030-07



Bracken Irrigation Cannula
 Mirror

Catalog # _____
 WI-030-17





Scalpel Handles
Skalpellgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Brierley LASIK Cannula
Mirror

Catalog # _____

WI-010-03



Buratto LASIK Cannula
Mirror

Catalog # _____

WI-010-04



Byrd Suction Tube
Satin

Catalog # _____

WI-010-34

Overall Length

21.5 cm – 8 1/2



Carter LASIK Cannula
Mirror

Catalog # _____

WI-012-03



Chevalier-Jackson Tracheal Tube
Mirror

Catalog # _____


WI-020-13



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Coakley Antrum Trocar
 Satin

Catalog #	Overall Length
WI-020-14	14 cm – 5 1/2



Cogswell Suction Tube
 Satin

Catalog #	Overall Length
WI-021-13	24.5 cm – 9 3/4



With Finger Cutt Off

Cohen Uterine Cannula
 Satin

Catalog #	Overall Length
WI-021-14	



Complete Cannula
 With 2 Cones

Cone Ventricular Cannula
 Satin

Catalog #	Overall Length
WI-030-07	



Cooley Suction Tube
 Satin

Catalog #	Overall Length
WI-030-17	36 cm – 14 1/4



36 cm, Ø 8.0 mm 31 cm, Ø 10.0 mm 31 cm, Ø 7.0 mm 36 cm, Ø 8.0 mm



Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Corbin Sub-Tenon Cannula
Mirror

Catalog # _____


WI-010-03



Coupland Suction Tube
Satin

Catalog # _____

WI-010-04



DeBaKey Suction Tube
Satin

Catalog # _____

WI-010-34



Drews Irrigating Vectis
Mirror

Catalog # _____

WI-012-03



Drof Infusion Cannula
Mirror

Catalog # _____

WI-020-13



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Eicken-Killian Irrigation Cannula
 Satin

Catalog #	Overall Length
WI-020-14	14 cm – 5 1/2



Fergusson Suction Tube
 Satin

Catalog #	Overall Length
WI-021-13	19.5 cm – 7 3/4



Finsterer Suction Tube
 Satin

Catalog #	Overall Length
WI-021-14	32 cm – 12 3/4



Frazier Suction Tube
 Satin

Catalog #	Overall Length
WI-030-07	19.5 cm – 7 3/4



Frazier Ventricular Cannula
 Satin

Catalog #	Overall Length
WI-030-17	9 cm – 3 1/2



Graduated in centimeters



Scalpel Handles
Skalpеллgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Fukasaku Anesthesia Cannula
Mirror

Catalog # _____

WI-010-03



Fukasaku Hydrochop Cannula
Mirror

Catalog # _____


WI-010-04



Fukasaku LASIK Cannula
Mirror

Catalog # _____

WI-010-34



Fukushima Tear Drop Suction Tube
Satin

Catalog # _____

WI-012-03



Gills Aspirating Cannula
Mirror

Catalog # _____

WI-020-13



Gills Irrigation Aspirating Cannula
Mirror

Catalog # _____
WI-020-14

Scalpel Handles
Skalpellgriffe
Mangos de bisturi
Manches pour bisturi
Manici per bisturi



Gimbel Cortex Aspirating Cannula
Mirror

Catalog # _____
WI-021-13



Gorney (Freer) Suction Elevator
Satin

Catalog #	Overall Length
WI-021-14	19 cm – 7 1/2



Gruber Suction Elevator
Satin

Catalog #	Overall Length
WI-030-07	14.5 cm – 5 3/4



Guilen Suction Raspatory
Satin

Catalog #	Overall Length
WI-030-17	19 cm – 7 1/2





Scalpel Handles
Skalpelli
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Haraldsson Suction Elevator
Satin

Catalog #	Overall Length
WI-010-03	18.5 cm – 7 1/4



Hartmann Ear Catheter
Satin

Catalog #	Overall Length
WI-010-04	15 cm – 6



Heslin Infusion Cannula
Mirror

Catalog #	Overall Length
WI-010-34	



House Suction Adapter
Satin

Catalog #	Overall Length
WI-012-03	



With Finger Cut Off

House Suction Tube
Satin

Catalog #	Overall Length
WI-020-13	



With Finger Cut Off



Scalpel Handles
 Skalpellgriffe
 Mangos de bisturíes
 Manches pour bistouris
 Manici per bisturi

Angled at 8 mm
 Flat Tip, 25 G

Hydrodissection Cannula
 Mirror

Catalog # _____
 WI-020-14



Infusion Cannula
 Mirror

Catalog # _____
 WI-021-13



Self Retaining Tip
 22 G

Irrigation Cannula
 Mirror

Catalog # _____
 WI-021-14



Irrigation Handpiece
 Mirror

Catalog # _____ Overall Length
 WI-030-07 10.5 cm – 4



Jensen Capsule Polisher Cannula
 Mirror

Catalog # _____
 WI-030-17





Scalpel Handles
Skalpgriffe
Mangos de bisturíes
Manches pour bistouris
Manici per bisturi

Amalgam Carrier
Satin

Catalog # _____

WI-010-03



Bone Graft Syringe
Satin

Catalog # _____

WI-010-04



Cook Waite Dental Aspirating Syringe
Satin

Catalog # _____

WI-010-34



Dental Aspirating Syringe
Satin

Catalog # _____

WI-012-03



Ear Syringe
Satin

Catalog # _____

WI-020-13





With Protection Disk and 2 Tips
50ml/100ml



Ivory Amalgam Carrier
Satin

Catalog #
WI-020-14






WAQAR INTERNATIONAL

Precision Instruments Manufacturer

Premium Quality Instruments Manufacturer Of Surgical, Dental, Veterinary Instruments.



OMAN OFFICE:

Phone: +968 9368 9750

Email: sales@waqarintl.com

Address: Muscat, Oman

HEAD OFFICE PAKISTAN

Phone: +92 307-6173627

Email: sales@waqarintl.com

Address: Roras Road, Opposite Anwar Khawaja Industries, Sialkot - 51310 Pakistan.

